

Profile Definition

MEDIBUS.X

WARNING

For a full understanding of this profile definition, the user should carefully read this document as well as the Instructions for Use of the basic device.

**For data communication
V 1.n**

Safety

This document contains the MEDIBUS response codes to be used by devices, which comply to the Common Medibus Device Profile MEDIBUS.X. Strictly follow the Instructions for Use of the corresponding therapy device using MEDIBUS.X profile. Any use of the MEDIBand US protocol and the MEDIBUS.X profile requires full understanding and strict observation of this document and the Instructions for Use of the mentioned medical devices, as appropriate. For a general description of the protocol, refer to Dräger RS 232 MEDIBUS Protocol Definition (order-no. 9028258, 13th edition 2016). For a general description of the MEDIBUS.X profile refer to "Dräger MEDIBUS.X Rules and Standards for Implementation" (order-no. 9052607, 5 – 2019-02).

Trademarks

- Evita
- Infinity
- Carina
- Babylog
- Savina
- Oxylog
- Atlan
- Zeus
- MEDIBUS.X
- Perseus
- Fabius
- Isolette
- Babyleo
- BiliLux
- Babyroo

are trademarks owned by Dräger.

Contents

Sheet	Description
Device identification numbers	List of devices supporting MEDIBUS.X together with their MEDIBUS ID Number.
Issued commands	List of commands, which are issued by Dräger devices.
Responded commands	List of commands, which are processed by Dräger devices. Commands not listed will be answered as specified for unexpected commands, but the device will not process the command any further.
Device setting response	List of response codes for the REQUEST CURRENT DEVICE SETTING command.
Current measured data and alarm limit response	List of response codes for the REQUEST CURRENT MEASURED DATA and REQUEST CURRENT ALARM LIMIT commands. The columns LL and HL indicate all devices that use the respective code for sending the upper or lower alarm limit.
Real-time data	List of supported real time codes.
Alarm status response	List of response codes for the REQUEST CURRENT ALARMS commands.
Text messages	List of response codes for the REQUEST CURRENT TEXT MESSAGES commands.
Infusion pump setting responses	List of response codes for the REQUEST CURRENT DEVICE SETTINGS sub-command of the DEVICE EXTENSION COMMAND FOR INFUSION PUMPS
Infusion pump current measured data and alarm limit response	List of response codes for the REQUEST MEASURED DATA and REQUEST ALARM LIMIT sub-command of the DEVICE EXTENSION COMMAND FOR INFUSION PUMPS. The columns LL and HL indicate all devices that use the respective code for sending the upper or lower alarm limit.
Infusion pump alarm status response	List of response codes for the REQUEST CURRENT ALARM sub-command of the DEVICE EXTENSION COMMAND FOR INFUSION PUMPS
Infusion pump text messages	List of response codes for the REQUEST CURRENT TEXT MESSAGES sub-command of the DEVICE EXTENSION COMMAND FOR INFUSION PUMPS

Document history

Edition	Description
Edition 22	This edition covers MEDIBUS.X 01.15. The new software 3.0 for Evita V600, Evita V800, Babylog VN600 and Babylog VN800 has been considered with this edition. See the code list documentation for further details.
Edition 21	This edition covers MEDIBUS.X 01.14. The new software 2.10 for Atlan has been considered with this edition. The former edition did not document that Atlan is no longer supporting a lower alarm limit for the measurement with code 05H on code page 1 in software 2.00. This has been corrected. The former edition did not document that Atlan is no longer supporting the alarm with code DAH on code page 1 in software 2.00. This has been corrected. The alarm code D3H on page 1 is used by Atlan software 1.00 and later. The code has not been documented in former editions of this document. This has been corrected. The new software 5.03 for Savina 300 has been considered with this edition. See the code list documentation for further details.
Edition 20	This edition covers MEDIBUS.X 01.13. The new software 2.00 for Atlan has been considered with this edition. The new software 2.00 for Evita V600, Evita V800, Babylog VN600 and Babylog VN800 has been considered with this edition. The new product Babyroo TN300 has been considered with this edition. Former editions falsely documented that Atlan is supporting text 5D as of software version 1.02. This has been corrected. See the code list documentation for further details.
Edition 19	This edition covers MEDIBUS.X 1.12. The new software 2.00 for Oxylog VE300 has been considered with this edition. The new software 2.02 for Zeus IE has been considered with this edition. The new software 2.03 for Perseus A500 has been considered with this edition. These software updates required to add response code definitions to the MEDIBUS.X profile. The new software 1.02 for Atlan has been considered with this edition. The new products Evita V600, Evita V800, Babylog VN600 and Babylog VN800 have been considered with this edition. These new devices use a revised nomenclature for some parameters, see the code list documentation for further details. Former editions did not document that ZEUS IE is supporting alarm E5 on code page 2 as of software version 01.10. This has been corrected.
Edition 18	This edition covers MEDIBUS.X 1.11. The new product Atlan has been considered with this edition. Clarification that gas consumption data is accumulated on a case by case basis. Some codes have incorrectly been associated with the Primus Family or its members. These associations have been removed in this edition. See the code list documentation for further details.

Edition	Description
Edition 17	The new software 1.10 for Oxylog VE300 has been considered with this edition. The version of the MEDIBUS.X profile remains 1.10.
Edition 16	The new software 2.01 for Zeus IE has been considered with this edition. The version of the MEDIBUS.X profile remains 1.10.
Edition 15	This edition covers MEDIBUS.X 1.10 and contains some minor editorial corrections. The new software 2.00 for Perseus A500 has been considered with this edition. The new products Oxylog VE300 and BiliLux have been considered with this edition. See the code list documentation for further details.
Edition 14	This edition covers MEDIBUS.X 1.09 The new software 1.00 for the Babyleo TN500 has been considered with this edition. See the code list documentation for further details.
Edition 13	This edition covers MEDIBUS.X 1.08 The new software 2.00 for the Zeus IE has been considered with this edition. The new software 5.00 for the Savina 300 has been considered with this edition. The new software 5.00 for the Isolette 8000 plus has been considered with this edition. This version is the first Isolette 8000 plus software version that supports MEDIBUS.X. See the code list documentation for further details.
Edition 12	This edition covers MEDIBUS.X 1.07 The new software 4.52 for the Primus family has been considered with this edition.
Edition 11	This edition covers MEDIBUS.X 1.06. The new software 2.42 for Evita Infinity V500, Evita V300 and Babylog VN500 has been considered with this edition. See the code list documentation for further details.
Edition 10	The Fabius Family has been added as a user of the MEDIBUS.X profile. The version of the profile remains 1.05.
Edition 9	This edition covers MEDIBUS.X 1.05. New codes have been added for the latest software releases of Zeus IE, Savina 300 and for future use. With this software release, the Savina also supports text messages for CO2 units. See the code list documentation for further details.
Edition 8	Not published.
Edition 7	This edition covers MEDIBUS.X 1.04. The code lists have been updated to support the latest software versions for MEDIBUS.X compliant devices. See the code list documentation for further details. The edition contains some further corrections of the former documentation for the Primus family. See the code list documentation for further details.
Edition 6	This edition covers MEDIBUS.X 1.03. The code lists have been extended to support the latest software versions for MEDIBUS.X compliant devices. See the code list documentation for further details. Other editorial changes: The document is now published as Excel file, minor editorial changes were necessary to support this new format. The commands 2C - Request current low alarm limits (codepage 2) and 2D Request current high alarm Limits (codepage 2) are supported by the Evita V500 and the Babylog VN500. This information was missing in the last edition. Same applies to the command 4A Configure data response command. The code lists now reflect that not all devices of the Primus family send out all codes defined for them. For other corrections see the information given in the following tables.

Edition
Edition 5

Description

This edition covers MEDIBUS.X 1.02.

The code lists have been extended to support the latest software versions for MEDIBUS.X compliant devices. See the code list documentation for further details.

Other editorial changes:

The documentation for limit codes gives information about the devices using the limit code.

The former edition did not document that the measure codes CP1 88H and CP2 22H are being used to send a lower limit value.

The text message 6CH is currently not send out by any device.

For Babylog VN500:

The former edition stated that Babylog VN500 transmits setting code 17H, measure codes CP1 0EH, 10H, 11H, CP2 25H and text message codes 0EH, 1DH, 5EH, but the device actually does not transmit these codes.

For Oxylog 3000 plus:

The former edition stated that Oxylog 3000 plus transmits alarm code CP1 42H, but the device actually does not transmit this code.

For Carina:

The former edition stated that Carina transmits measure code CP1 B5, but the device actually does not transmit this code.

The former edition stated that Carina transmits command real-time configuration changed 56H, but the device actually does not transmit this command.

The former edition did not document that Carina transmits the measure code CP1 06H and D6H, CP2 22H and the setting codes 42H and 44H.

Carina SW 3.20 also transmits codes that have been introduced with profile version 1.02. Details see below.

Edition 4

Code mapping for Perseus A500 added.

Alarm code 3E (page 1) is not used by Savina 300

Edition 3

Supports the following devices: Zeus IE, Savina 300, Primus family, Oxylog 3000 plus, Carina, Evita V500, Babylog VN500

Code list extended to support these devices and the latest software versions for other MEDIBUS.X compliant devices.

Section for infusion pump support added.

Editorial changes and spelling corrections for existing code descriptions.

Adaptions in parameter nomenclature.

Deprecate use of some reserved codes that were introduced as a preparatory step for legacy reasons with 01.00.

Edition 2

Description for measure code 04 (code page 1) emended.

Device identification numbers

Device	ID number	First SW version that supports Medibus.X	Remarks
Infinity Acute Care System Evita Infinity V500	5000	SW 2.20	
Infinity Acute Care System Babylog VN500	5001	SW 2.20	
Evita V300	5002	SW 2.31	Edition 6: New device Evita V300 supports MEDIBUS.X.
Evita V600	5004	SW 1.05	Edition 19: Evita V600 supports MEDIBUS.X beginning with SW 1.05
Babylog VN600	5005	SW 1.05	Edition 19: Babylog VN600 supports MEDIBUS.X beginning with SW 1.05
Evita V800	5006	SW 1.05	Edition 19: Evita V800 supports MEDIBUS.X beginning with SW 1.05
Babylog VN800	5007	SW 1.05	Edition 19: Babylog VN800 supports MEDIBUS.X beginning with SW 1.05
Savina 300	5030	SW 4.00	
Carina	5040	SW 3.20	
Oxylog 3000 plus	5050	SW 1.02	
Oxylog VE300	5051	SW 1.00	Edition 15: Oxylog VE300 supports MEDIBUS.X beginning with SW 1.00
Perseus A500	5500	SW 1.00	
Zeus Infinity Empowered	5501	SW 1.01	
Primus	5510	SW 4.50	
Primus Infinity Empowered	5511	SW 4.50	
Apollo	5512	SW 4.50	
Atlan	5555	SW 1.00	Edition 18: Atlan supports MEDIBUS.X beginning with SW 1.00
Fabius GS	5599	SW 3.35	Edition 10: The Fabius family supports MEDIBUS.X beginning with SW 3.35.
Fabius GS Premium	5600	SW 3.35	
Fabius OS	5601	SW 3.35	
Fabius Tiro, Fabius D-M	5602	SW 3.35	
Fabius MRI	5603	SW 3.35	
Fabius plus	5604	SW 3.35	
Fabius Plus XL	5605	SW 3.35	
Babyleo TN500	5700	SW 1.00	Edition 14: Babyleo TN500 supports MEDIBUS.X beginning with SW 1.00
Isolette 8000 plus	5701	SW 5.00	Edition 13: Isolette 8000 plus supports MEDIBUS.X beginning with SW 5.00

Device	ID number	First SW version that supports Medibus.X	Remarks
Babyroo TN300	5702	SW 1.00	Edition 20: Babyroo TN300 supports MEDIBUS.X beginning with SW 1.00
BiliLux	6300	SW 1.00	Edition 15: BiliLux supports MEDIBUS.X beginning with SW 1.00
Note: In the following lists, the term Fabius Family refers to all Fabius variants.			
The term Primus Family refers to Primus, Primus Infinity Empowered, and Apollo.			
The term Atlan refers to all product variants.			

Issued commands

Code	Code description	Used by device	Remarks
24H	Request current measured data (codepage 1)	Fabius Family	
28H	Request current date and time	Perseus A500, Atlan	
29H	Request current device settings	Perseus A500, Atlan	
2AH	Request current text messages	Perseus A500, Atlan, Fabius Family	
30H	No Operation (NOP)	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300, Carina, Babyleo TN500, Isolette 8000 plus, Babyroo TN300, Fabius Family, BiliLux	

Code	Code description	Used by device	Remarks
49H	Time changed	Perseus A500, Atlan, Evita V600, Evita V800, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300, Carina, Babyleo TN500, Babyroo TN300 (added with SW 3.00 for Evita V800, Evita V600, Babylog VN800 and Babylog VN600)	Edition 22: Evita V600, Evita V800, Babylog VN600 and Babylog VN800 issue this command as of SW 3.00.
4AH	Configure data response command	Perseus A500, Atlan	
51H	Initialize Communication (ICC)	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300, Carina, Babyleo TN500, Isolette 8000 plus, Babyroo TN300, Fabius Family, BiliLux	

Code	Code description	Used by device	Remarks
52H	Request Device Identification	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300, Carina, Babyleo TN500, Isolette 8000 plus, Babyroo TN300, Fabius Family, BiliLux	
55H	Stop Communication (STOP)	Perseus A500, Atlan	
56H	Real-time configuration changed	Perseus A500, Atlan, Primus Family, Zeus IE	

Responded commands

Code	Code description	Used by device	Remarks
24H	Request current measured data (codepage 1)	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300, Carina, Babyleo TN500, Isolette 8000 plus, Babyroo TN300, Fabius Family, BiliLux	
25H	Request current low alarm limits (codepage 1)	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300, Carina, Fabius Family	

Code	Code description	Used by device	Remarks
26H	Request current high alarm Limits (codepage 1)	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300, Carina, Fabius Family	
27H	Request current alarms (codepage 1)	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300, Carina, Babyleo TN500, Isolette 8000 plus, Babyroo TN300, Fabius Family	

Code	Code description	Used by device	Remarks
28H	Request current date and time	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300, Carina, Babyleo TN500, Babyroo TN300, Fabius Family	
29H	Request current device settings	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300, Carina, Babyleo TN500, Isolette 8000 plus, Babyroo TN300, Fabius Family	

Code	Code description	Used by device	Remarks
2AH	Request current text messages	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300, Carina, Babyleo TN500, Isolette 8000 plus, Babyroo TN300, Fabius Family, BiliLux	
2BH	Request current measured data (codepage 2)	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Carina, Babyleo TN500	
2CH	Request current low alarm limits (codepage 2)	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Zeus IE	
2DH	Request current high alarm Limits (codepage 2)	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Zeus IE	

Code	Code description	Used by device	Remarks
2EH	Request current alarms (codepage 2)	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300, Carina, Babyleo TN500, Isolette 8000 plus, Babyroo TN300, Fabius Family	
30H	No Operation (NOP)	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300, Carina, Babyleo TN500, Isolette 8000 plus, Babyroo TN300, Fabius Family, BiliLux	
49H	Time changed	Perseus A500, Atlan	

Code	Code description	Used by device	Remarks
4AH	Configure data response command	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300, Carina, Babyleo TN500, Isolette 8000 plus, Babyroo TN300	
51H	Initialize Communication (ICC)	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300, Carina, Babyleo TN500, Isolette 8000 plus, Babyroo TN300, Fabius Family, BiliLux	

Code	Code description	Used by device	Remarks
52H	Request Device Identification	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300, Carina, Babyleo TN500, Isolette 8000 plus, Babyroo TN300, Fabius Family, BiliLux	
53H	Request real-time configuration	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300, Carina, Babyroo TN300, Fabius Family	

Code	Code description	Used by device	Remarks
54H	Configure real-time transmission	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300, Carina, Babyroo TN300, Fabius Family	
55H	Stop Communication (STOP)	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300, Carina, Babyleo TN500, Isolette 8000 plus, Babyroo TN300, Fabius Family, BiliLux	
59H	Device extension command for infusion pumps	Zeus IE	

Device setting response

Code	Data description	Unit	New with version	Used by device	Remarks
01H	Inspiratory oxygen fraction, FiO2	%	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus, Primus Infinity Empowered, Zeus IE, Oxylog 3000 plus, Carina (removed with SW 4.51 for Primus, Primus Infinity Empowered)	Edition 7: Members of the Primus family used codes in their legacy meaning. This has been corrected with SW 4.51
02H	Inspiratory flow, Flowi	L/min	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Fabius Family	
03H	Time interval between sighs, Interval sigh	s	01.01	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800	

Code	Data description	Unit	New with version	Used by device	Remarks
04H	Inspiratory tidal volume, VTi	L	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300, Carina, Fabius Family	
05H	Inspiratory time, Ti	s	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300, Carina, Fabius Family	
06H	Expiratory time, Te	s	01.01	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Zeus IE	

Code	Data description	Unit	New with version	Used by device	Remarks
07H	I : E I part	no unit	01.00	Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Zeus IE, Fabius Family (added with SW 2.00 for Atlan)	
08H	I : E E part	no unit	01.00	Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Zeus IE, Fabius Family (added with SW 2.00 for Atlan)	
09H	Respiratory rate, RR	1/min	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300, Carina, Fabius Family	

Code	Data description	Unit	New with version	Used by device	Remarks
0AH	Minimal provided respiratory rate, RRmin	1/min	01.05	Perseus A500, Atlan, Primus Family, Zeus IE, Fabius Family (added with SW 2.00 for Perseus A500, added with SW 1.10 for Zeus IE, added with SW 4.52 for the Primus Family)	
0BH	Positive end-expiratory pressure, PEEP	mbar	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300, Carina, Fabius Family (added with SW 4.51 for the Primus Family)	Edition 7: Members of the Primus family used codes in their legacy meaning. This has been corrected with SW 4.51

Code	Data description	Unit	New with version	Used by device	Remarks
0CH	Additional intermittent PEEP for sighs, $\Delta \text{intPEEPsigh}$	mbar	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE (removed with SW 1.10 for Zeus IE, removed with SW 4.51 for the Primus Family)	Edition 9: Zeus IE used this code in its legacy meaning. This has been corrected with SW 1.10 Edition 7: Members of the Primus family used codes in their legacy meaning. This has been corrected with SW 4.51
0DH	Lower pressure level during APRV, P _{low}	mbar	01.00	Perseus A500, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300	
0EH	Upper pressure level during APRV, P _{high}	mbar	01.00	Perseus A500, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300	
0FH	Time of lower pressure level during APRV, T _{low}	s	01.00	Perseus A500, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300	

Code	Data description	Unit	New with version	Used by device	Remarks
10H	Time of upper pressure level during APRV, Thigh	s	01.00	Perseus A500, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300	
11H	Apnea alarm time, Tapn	s	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Oxylog 3000 plus, Oxylog VE300, Carina (added with SW 2.00 for Oxylog VE300)	
12H	Pressure amplitude above PEEP in Pressure Support, ΔP_{supp}	mbar	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300, Carina, Fabius Family	

Code	Data description	Unit	New with version	Used by device	Remarks
13H	Pressure limitation, Pmax	mbar	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300, Carina, Fabius Family	
14H	Number of respiratory cycles per sigh phase, Cycles sigh	cycles/sigh	01.01	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800	Edition 19: Clarification of parameter semantics.
17H	P01 repetition interval, Interval P0.1	min	01.01	Evita V300, Evita Infinity V500, Evita V600, Evita V800	
18H	Delay for minute volume high alarm, MVhigh delay	s	01.01	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800	
19H	Delay for minute volume low alarm, MVlow delay	s	01.01	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800	
1AH	Air temperature, T Air	°C	01.08	Babyleo TN500, Isolette 8000 plus	

Code	Data description	Unit	New with version	Used by device	Remarks
1BH	Skin temperature, T Skin	°C	01.08	Babyleo TN500, Isolette 8000 plus, Babyroo TN300	
1CH	Relative humidity, Humidity	%	01.08	Babyleo TN500, Isolette 8000 plus	
24H	Oxygen concentration, O2	%	01.08	Babyleo TN500, Isolette 8000 plus	
27H	Plateau time in % of inspiratory time, % Tplat	%	01.00	Perseus A500, Atlan, Primus Family, Zeus IE, Oxylog 3000 plus, Fabius Family	
29H	Flow trigger threshold, Trigger (Synonymous parameter label: Flow trigger)	L/min	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300, Fabius Family	Edition 19: New devices use a different terminology for the same parameter.
2AH	Frequency during HFO, fhf	Hz	01.00	Babylog VN500, Babylog VN600, Babylog VN800	
2EH	Pressure rise time, Slope	s	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300, Carina	

Code	Data description	Unit	New with version	Used by device	Remarks
2FH	Fresh-gas flow, FG flow	mL/min	01.00	Perseus A500, Atlan, Primus, Primus Infinity Empowered, Zeus IE	Edition 6: Apollo does not send this code.
3CH	Flow assist during PPS ventilation, Flow Assist	(mbar*s/L)/10	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800	Edition 12: Introduced brackets to underline operator priority.
3DH	Volume assist during PPS ventilation, Vol. Assist	(mbar/L)/10	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800	Edition 12: Introduced brackets to underline operator priority.
3EH	Radiant warmer power, Heater rad	%	01.09	Babyleo TN500, Babyroo TN300	
3FH	Mattress temperature, T Matt	°C	01.09	Babyleo TN500, Babyroo TN300	
40H	Inspiratory tidal volume, VTi	mL	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Zeus IE	

Code	Data description	Unit	New with version	Used by device	Remarks
42H	Respiratory rate for apnea ventilation, RRapn	1/min	01.00	Perseus A500, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Oxylog 3000 plus, Oxylog VE300, Carina (removed with SW 2.00 for Perseus A500, removed with SW 4.52 for the Primus Family, added with SW 2.00 for Oxylog VE300)	Edition 15: Perseus A500 used this code for the device parameter RRmin. This has been corrected with SW 2.00 Edition 12: Members of the Primus family used this code for the device parameter RRmin. This has been corrected with SW 4.52
44H	Tidal volume for apnea ventilation, VTapn	L	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Oxylog 3000 plus, Oxylog VE300, Carina (added with SW 2.00 for Oxylog VE300)	

Code	Data description	Unit	New with version	Used by device	Remarks
45H	Inspiratory pressure, P _{insp}	mbar	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300, Carina, Fabius Family (added with SW 1.10 for Oxylog VE300)	
46H	Compensation rate during ATC, % comp	%	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800	
47H	Inner diameter of the endotracheal/tracheostoma tube, Ø tube	mm	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300 (added with SW 5.00 for Savina 300)	
49H	Ratio of inspiration to expiration I:E	%	01.01	Zeus IE, Oxylog 3000 plus, Oxylog VE300, Carina	
4AH	Age	a	01.00	Perseus A500, Atlan, Primus Family, Zeus IE	

Code	Data description	Unit	New with version	Used by device	Remarks
4BH	Weight, Wt	kg	01.00	Perseus A500, Atlan, Zeus IE	
4EH	Disconnection alarm time, Tdiscon	s	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Carina	
4FH	Flow Acceleration, Flow Acc	mbar/s	01.00	Savina 300	
55H	Volatile anesthetic agent concentration in the fresh gas, FG VA	kPa	01.01	Zeus IE	
56H	End-tidal volatile agent concentration, etVA	kPa	01.01	Zeus IE	
5BH	Isoflurane concentration in the fresh gas, FG Iso	kPa	01.00	Perseus A500, Zeus IE	
5CH	End-tidal isoflurane concentration, etIso	kPa	01.01	Zeus IE	
5DH	Desflurane concentration in the fresh gas, FG Des	kPa	01.00	Perseus A500, Zeus IE	
5EH	End-tidal desflurane concentration, etDes	kPa	01.01	Zeus IE	
5FH	Sevoflurane concentration in the fresh gas, FG Sev	kPa	01.00	Perseus A500, Zeus IE	
60H	End-tidal sevoflurane concentration, etSev	kPa	01.01	Zeus IE	
61H	Oxygen concentration in the fresh gas, FG O2	%	01.00	Perseus A500, Atlan, Primus, Primus Infinity Empowered, Zeus IE (added with SW 4.51 for Primus, Primus Infinity Empowered)	Edition 7: Members of the Primus family used codes in their legacy meaning. This has been corrected with SW 4.51
62H	Minimum fresh-gas flow, MinFG flow	L/min	01.01	Zeus IE	
63H	Maximal inspiratory time for flow cycled strokes, Timax	s	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Carina	

Code	Data description	Unit	New with version	Used by device	Remarks
65H	Constant flow during O2 therapy, Const flow (Synonymous parameter label: Flow)	L/min	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300 (added with SW 5.00 for Savina 300)	Edition 19: Clarification of parameter semantics. New devices use a different terminology for the same parameter.
66H	Mean airway pressure during HFO, MAPhf	mbar	01.00	Babylog VN500, Babylog VN600, Babylog VN800	
67H	Inspiratory part of I:E ratio during HFO, I (I:Ehf)	no unit	01.00	Babylog VN500, Babylog VN600, Babylog VN800	Edition 19: Clarification of parameter semantics.
68H	Expiratory part of I:E ratio during HFO, E (I:Ehf)	no unit	01.00	Babylog VN500, Babylog VN600, Babylog VN800	Edition 19: Clarification of parameter semantics.
69H	Pressure amplitude during HFO, Ampl hf	mbar	01.00	Babylog VN500, Babylog VN600, Babylog VN800	
6AH	Tidal volume during HFO, VThf	mL	01.00	Babylog VN500, Babylog VN600, Babylog VN800	
6BH	Respiratory rate of sighs during HFO, RR _{sigh}	1/min	01.00	Babylog VN500, Babylog VN600, Babylog VN800	
6CH	Inspiratory time of sighs during HFO, Tisigh	s	01.00	Babylog VN500, Babylog VN600, Babylog VN800	
6DH	Inspiratory pressure of sighs during HFO, Psigh	mbar	01.00	Babylog VN500, Babylog VN600, Babylog VN800	
6EH	Tidal volume during apnea ventilation, VT _{apn}	mL	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800	
6FH	Ideal body weight, IBW	kg	01.01	Primus Family, Zeus IE (added with SW 2.01 for Zeus IE)	

Code	Data description	Unit	New with version	Used by device	Remarks
72H	Inspiratory termination criterion based on peak inspiratory flow, Insp. term.	%PIF	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Savina 300, Oxylog 3000 plus	
73H	Expiratory termination criterion during APRV/AutoRelease, Exp. term.	%PEF	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800	
74H	Ideal body weight in SmartCare, SC-Weight (synonymous parameter label: SC-Wt)	kg	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800	Edition 19: New devices use a different terminology for the same parameter.
75H	Start of night rest in SmartCare, SC-start of rest (synonymous parameter label: SC-startrest)	h	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800	Edition 19: New devices use a different terminology for the same parameter.
76H	End of night rest in SmartCare, SC-end of rest (synonymous parameter label: SC-endrest)	h	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800	Edition 19: New devices use a different terminology for the same parameter.
78H	Pressure rise time of sighs during HFO, Slopesigh	s	01.02	Babylog VN500, Babylog VN600, Babylog VN800 (added with SW 2.30 for Babylog VN500)	Edition 19: Clarification of parameter semantics.
79H	FiO2 in SmartCare - upper limit, SC-FiO2 high	%	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800	Edition 19: Parameter label correction
7AH	PEEP in SmartCare - upper limit, SC-PEEPmax	mbar	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800	
7BH	Pressure variability, Press. var.	%	01.00	Evita Infinity V500, Evita V600, Evita V800	

Code	Data description	Unit	New with version	Used by device	Remarks
7DH	Maximum time of low pressure level during APRV/AutoRelease, Tlow max	s	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800	
7EH	Spontaneous respiratory rate in SmartCare - lower limit, SC-RRspon low	1/min	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800	
7FH	Spontaneous respiratory rate in SmartCare - upper limit, SC-RRspon high	1/min	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800	
80H	Tidal volume in SmartCare - lower limit, SC-VT low	mL/kg	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800	Edition 19: Parameter label correction
81H	End-tidal CO2 in SmartCare - upper limit, SC-etCO2high	mmHg	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800	
82H	Patient height in SmartCare, SC-Height (synonymous parameter label: SC-Ht)	cm	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800	Edition 19: New devices use a different terminology for the same parameter.
83H	Patient height, Ht	cm	01.00	Perseus A500, Atlan, Zeus IE (added with SW 2.01 for Zeus IE)	
84H	Maximal peak to peak pressure during HFO, Ampl hf max	mbar	01.00	Babylog VN500, Babylog VN600, Babylog VN800	
85H	Volume limitation, VTmax	mL	01.00	Babylog VN500, Babylog VN600, Babylog VN800 (removed with SW 3.00 for Babylog VN800 and Babylog VN600)	

Code	Data description	Unit	New with version	Used by device	Remarks
86H	Manual inspiratory pressure for Neonates during NIV - Non-Invasive Ventilation, PmanInsp	mbar	01.02	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800 (added with SW 2.30 for Babylog VN500, added with SW 2.30 for Evita Infinity V500)	
87H	Maximum time of pressure for manual inspiration pressure for Neonates during NIV - Non-Invasive Ventilation, TmanInsp	s	01.02	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800 (added with SW 2.30 for Babylog VN500, added with SW 2.30 for Evita Infinity V500)	
88H	Maximum airway pressure in Smart Ventilation Control, SVC-Pmax	mbar	01.08	Zeus IE (added with SW 2.00 for Zeus IE)	
8AH	Airway access Enumeration value with the following definition: 0: Tube, 1: Laryngeal Mask	no Unit	01.08	Zeus IE (added with SW 2.00 for Zeus IE, removed with SW 2.01 for Zeus IE)	
8BH	Lung mechanics in Smart Ventilation Control, Lung mechanics Enumeration value with the following definition: 0: Normal, 1: Restrictive, 2: Obstructive	no Unit	01.08	Zeus IE (added with SW 2.00 for Zeus IE)	
8CH	Ventilation goal of Smart Ventilation Control, Ventilation goal Enumeration value with the following definition: 0: Controlled Ventilation, 1: Augmented Ventilation, 2: Encourage SB, 3: Recovery	no Unit	01.08	Zeus IE (added with SW 2.00 for Zeus IE)	

Current measured data and alarm limit response

Code	Code page	Data description	Unit	LL	HL	New with version	Used by device	Remarks
00H	1	Accumulated air consumption of the device per case, Air cons (Note: 9999 L equals to 9999 L or more.)	L			01.00	Perseus A500, Atlan, Zeus IE	Edition 18: Clarification: accumulation happens per case, not per device lifetime.
01H	1	Accumulated N2O consumption of the device per case, N2O cons (Note: 9999 L equals to 9999 L or more.)	L			01.00	Perseus A500, Atlan, Zeus IE	Edition 18: Clarification: accumulation happens per case, not per device lifetime.
03H	1	Expiratory minute volume, MVe	L/min	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300	01.01	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300 (added with SW 5.03 for Savina 300)	
04H	1	Accumulated O2 consumption of the device per case, O2 cons (Note: 9999 L equals to 9999 L or more.)	L			01.00	Perseus A500, Atlan, Zeus IE	Edition 18: Clarification: accumulation happens per case, not per device lifetime.
05H	1	Airway pressure, Paw (Code is not used as response to GET MEASURED DATA command, only for sending alarm limits)	mbar	Perseus A500, Primus Family	Perseus A500, Atlan, Primus Family	01.00	Perseus A500, Atlan, Primus Family	Edition 21: As of SW 2.00 Atlan does not send this code within responses to the command 25H "Request current low Alarm Limits (Codepage 1)"
06H	1	Inspiratory minute volume, MVi	L/min	Carina	Carina	01.01	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Carina	

Code	Code page	Data description	Unit	LL	HL	New with version	Used by device	Remarks
07H	1	Dynamic compliance, Cdyn	L/bar			01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE	
08H	1	Resistance, R	mbar/(L/s)			01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Zeus IE	Edition 15: Specification of unit now contains brackets.
09H	1	CO2 production, V'CO2	mL/min			01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN600, Babylog VN800, Primus Family, Zeus IE (added with SW 2.00 for Zeus IE)	

Code	Code page	Data description	Unit	LL	HL	New with version	Used by device	Remarks
0AH	1	Spontaneous expiratory minute volume, MVespon	L/min			01.01	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family (added with SW 5.03 for Savina 300)	
0BH	1	Patient airway resistance, Rpat	mbar/(L/s)			01.01	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800	Edition 15: Specification of unit now contains brackets.
0DH	1	Mandatory expiratory minute volume, MVemand	L/min			01.01	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Primus Family	
0EH	1	Increase of measured CO2 value in phase III of the capnogram, CO2 slope	mmHg/L			01.01	Evita V300, Evita Infinity V500, Evita V600, Evita V800	.
0FH	1	O2 central gas supply pressure, O2-CGS	mbar			01.01	Zeus IE	
10H	1	Increase of measured CO2 value in phase III of the capnogram, CO2 slope	kPa/L			01.01	Evita V300, Evita Infinity V500, Evita V600, Evita V800	
11H	1	Increase of measured CO2 value in phase III of the capnogram, CO2 slope	Vol%/L			01.01	Evita V300, Evita Infinity V500, Evita V600, Evita V800	

Code	Code page	Data description	Unit	LL	HL	New with version	Used by device	Remarks
12H	1	Average device flow, Flow device (synonymous parameter label: FlowDev)	L/min			01.01	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800	Edition 19: New devices use a different terminology for the same parameter.
16H	1	Mandatory tidal volume, VTmand	mL			01.02	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800 (added with SW 2.30 for Babylog VN500, added with SW 2.30 for Evita Infinity V500)	
19H	1	ST deviation in lead I, STI This value is transmitted without decimals and may be rounded compared to the displayed value. See also code 29H on code page 2.	mm			01.01	Zeus IE	Edition 9: Zeus IE is sending this parameters with two formats. The profile specifies the precision of this value to explain the difference.
1AH	1	ST deviation in lead II, STII This value is transmitted without decimals and may be rounded compared to the displayed value. See also code 2AH on code page 2.	mm			01.01	Zeus IE	Edition 9: Zeus IE is sending this parameters with two formats. The profile specifies the precision of this value to explain the difference.
1BH	1	Accumulated halothane consumption of the device, Hal cons	mL			01.00	Perseus A500, Atlan, Primus Family	
1CH	1	Accumulated enflurane consumption of the device, Enf cons	mL			01.00	Perseus A500, Atlan, Primus Family	
1DH	1	Accumulated isoflurane consumption of the device, Iso cons	mL			01.00	Perseus A500, Atlan, Primus Family, Zeus IE	

Code	Code page	Data description	Unit	LL	HL	New with version	Used by device	Remarks
1EH	1	Accumulated desflurane consumption of the device, Des cons	mL			01.00	Perseus A500, Atlan, Primus Family, Zeus IE	
1FH	1	Accumulated sevoflurane consumption of the device, Sev cons	mL			01.00	Perseus A500, Atlan, Primus Family, Zeus IE	
21H	1	Cardiac index, CI	L/min/m ²			01.01	Zeus IE	
22H	1	Stroke volume, SV	mL			01.01	Zeus IE	
23H	1	Stroke index, SVI	mL/m ²			01.01	Zeus IE	
24H	1	Left ventricular stroke work index, LVSWI	g*m/m ²			01.01	Zeus IE	
25H	1	Right ventricular stroke work index, RVSWI	g*m/m ²			01.01	Zeus IE	
26H	1	Systemic vascular resistance, SVR	dyn*s/cm ⁵			01.01	Zeus IE	Edition 12: The unit has been corrected.
27H	1	Accumulated gas consumption of the device, Gas cons	L			01.01	Zeus IE	
28H	1	Therapy case duration, Tcase	min			01.00	Perseus A500, Atlan	
2AH	1	Correlation factor, r ²	no unit			01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800	
2BH	1	Left ventricular pressure - diastolic, LV D	mmHg	Zeus IE	Zeus IE	01.01	Zeus IE	
2CH	1	Left ventricular pressure - systolic, LV S	mmHg	Zeus IE	Zeus IE	01.01	Zeus IE	
2DH	1	Left ventricular pressure - mean, LV M	mmHg	Zeus IE	Zeus IE	01.01	Zeus IE	
2EH	1	Right ventricular pressure - diastolic, RV D	mmHg	Zeus IE	Zeus IE	01.01	Zeus IE	
2FH	1	Right ventricular pressure - systolic, RV S	mmHg	Zeus IE	Zeus IE	01.01	Zeus IE	
30H	1	Right ventricular pressure - mean, RV M	mmHg	Zeus IE	Zeus IE	01.01	Zeus IE	
31H	1	General diastolic pressure (channel 1), GP1 D	mmHg	Zeus IE	Zeus IE	01.01	Zeus IE	
32H	1	General systolic pressure (channel 1), GP1 S	mmHg	Zeus IE	Zeus IE	01.01	Zeus IE	
33H	1	General mean pressure (channel 1), GP1 M	mmHg	Zeus IE	Zeus IE	01.01	Zeus IE	
34H	1	General diastolic pressure (channel 2), GP2 D	mmHg	Zeus IE	Zeus IE	01.01	Zeus IE	
35H	1	General systolic pressure (channel 2), GP2 S	mmHg	Zeus IE	Zeus IE	01.01	Zeus IE	
36H	1	General mean pressure (channel 2), GP2 M	mmHg	Zeus IE	Zeus IE	01.01	Zeus IE	
37H	1	Intracranial pressure, ICP	mmHg	Zeus IE	Zeus IE	01.01	Zeus IE	
38H	1	Cerebral perfusion pressure, CPP	mmHg	Zeus IE	Zeus IE	01.01	Zeus IE	

Code	Code page	Data description	Unit	LL	HL	New with version	Used by device	Remarks
39H	1	Left atrial pressure, LA	mmHg	Zeus IE	Zeus IE	01.01	Zeus IE	
3CH	1	Mean dynamic compliance, Cdyn mean	L/bar			01.10	Perseus A500, Atlan (added with SW 2.00 for Perseus A500)	
50H	1	Inspiratory halothane concentration, inHal	kPa	Perseus A500, Atlan, Primus Family, Zeus IE	Perseus A500, Atlan, Primus Family, Zeus IE	01.00	Perseus A500, Atlan, Primus Family, Zeus IE	
51H	1	End-tidal halothane concentration, etHal	kPa			01.00	Perseus A500, Atlan, Primus Family, Zeus IE	
52H	1	Inspiratory enflurane concentration, inEnf	kPa	Perseus A500, Atlan, Primus Family, Zeus IE	Perseus A500, Atlan, Primus Family, Zeus IE	01.00	Perseus A500, Atlan, Primus Family, Zeus IE	
53H	1	End-tidal enflurane concentration, etEnf	kPa			01.00	Perseus A500, Atlan, Primus Family, Zeus IE	
54H	1	Inspiratory isoflurane concentration, inIso	kPa	Perseus A500, Atlan, Primus Family, Zeus IE	Perseus A500, Atlan, Primus Family, Zeus IE	01.00	Perseus A500, Atlan, Primus Family, Zeus IE	
55H	1	End-tidal isoflurane concentration, etIso	kPa			01.00	Perseus A500, Atlan, Primus Family, Zeus IE	
56H	1	Inspiratory desflurane concentration, inDes	kPa	Perseus A500, Atlan, Primus Family, Zeus IE	Perseus A500, Atlan, Primus Family, Zeus IE	01.00	Perseus A500, Atlan, Primus Family, Zeus IE	
57H	1	End-tidal desflurane concentration, etDes	kPa			01.00	Perseus A500, Atlan, Primus Family, Zeus IE	
58H	1	Inspiratory sevoflurane concentration, inSev	kPa	Perseus A500, Atlan, Primus Family, Zeus IE	Perseus A500, Atlan, Primus Family, Zeus IE	01.00	Perseus A500, Atlan, Primus Family, Zeus IE	
59H	1	End-tidal sevoflurane concentration, etSev	kPa			01.00	Perseus A500, Atlan, Primus Family, Zeus IE	

Code	Code page	Data description	Unit	LL	HL	New with version	Used by device	Remarks
5AH	1	Inspiratory concentration of primary agent	kPa			01.00	Perseus A500, Atlan, Primus Family, Zeus IE	
5BH	1	End-tidal concentration of primary agent	kPa			01.00	Perseus A500, Atlan, Primus Family, Zeus IE	
5CH	1	Inspiratory concentration of secondary agent	kPa			01.00	Perseus A500, Atlan, Primus Family, Zeus IE	
5DH	1	End-tidal concentration of secondary agent	kPa			01.00	Perseus A500, Atlan, Primus Family, Zeus IE	
5EH	1	ST deviation in lead III, STIII This value is transmitted without decimals and may be rounded compared to the displayed value. See also code 2BH on code page 2.	mm			01.01	Zeus IE	Edition 9: Zeus IE is sending this parameters with two formats. The profile specifies the precision of this value to explain the difference.
5FH	1	ST deviation in lead aVL, STaVL This value is transmitted without decimals and may be rounded compared to the displayed value. See also code 2CH on code page 2.	mm			01.01	Zeus IE	Edition 9: Zeus IE is sending this parameters with two formats. The profile specifies the precision of this value to explain the difference.
60H	1	ST deviation in lead aVF, STaVF This value is transmitted without decimals and may be rounded compared to the displayed value. See also code 2DH on code page 2.	mm			01.01	Zeus IE	Edition 9: Zeus IE is sending this parameters with two formats. The profile specifies the precision of this value to explain the difference.
61H	1	ST deviation in lead aVR, STaVR This value is transmitted without decimals and may be rounded compared to the displayed value. See also code 2EH on code page 2.	mm			01.01	Zeus IE	Edition 9: Zeus IE is sending this parameters with two formats. The profile specifies the precision of this value to explain the difference.

Code	Code page	Data description	Unit	LL	HL	New with version	Used by device	Remarks
62H	1	ST deviation in lead V, STV This value is transmitted without decimals and may be rounded compared to the displayed value. See also code 2FH on code page 2.	mm			01.01	Zeus IE	Edition 9: Zeus IE is sending this parameters with two formats. The profile specifies the precision of this value to explain the difference.
63H	1	ST deviation in lead V+, STV+ This value is transmitted without decimals and may be rounded compared to the displayed value. See also code 30H on code page 2.	mm			01.01	Zeus IE	Edition 9: Zeus IE is sending this parameters with two formats. The profile specifies the precision of this value to explain the difference.
64H	1	O2 uptake, V'O2	10 mL/min			01.00	Perseus A500, Atlan, Primus Family	
65H	1	Peripheral skin temperature, T Peri	°C			01.08	Babyleo TN500, Isolette 8000 plus, Babyroo TN300	
66H	1	Central skin temperature, T Centr	°C			01.08	Babyleo TN500, Isolette 8000 plus, Babyroo TN300	
68H	1	Spontaneous inspiratory tidal volume, VTispon	mL			01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800	
69H	1	Plateau time, Tplat	s			01.00	Savina 300	
6AH	1	Mattress temperature, T Matt	°C			01.09	Babyleo TN500, Babyroo TN300	
6BH	1	Ambient pressure, Pamb	mbar			01.00	Perseus A500, Atlan, Primus Family, Zeus IE	
6CH	1	Relative humidity, Humidity	%			01.08	Babyleo TN500, Isolette 8000 plus	
6DH	1	Air temperature, T Air	°C			01.08	Babyleo TN500, Isolette 8000 plus	
6EH	1	Radiant warmer power, Heater rad	%			01.09	Babyleo TN500, Babyroo TN300	

Code	Code page	Data description	Unit	LL	HL	New with version	Used by device	Remarks
6FH	1	Spontaneous inspiratory time, Tispon	s			01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800	
71H	1	Minimum pressure, Pmin	mbar			01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800	
72H	1	Occlusion pressure in P0.1 maneuver, P0.1	mbar			01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800	
73H	1	Mean airway pressure, Pmean	mbar			01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Carina, Fabius Family	
74H	1	Plateau pressure, Pplat	mbar			01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Fabius Family	
76H	1	Inspiratory peak flow, Flowipeak	mL/s			01.00	Savina 300	

Code	Code page	Data description	Unit	LL	HL	New with version	Used by device	Remarks
78H	1	Positive end-expiratory pressure, PEEP	mbar			01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Carina, Fabius Family	
79H	1	Intrinsic PEEP, PEEPi	mbar			01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800	
7BH	1	Mandatory respiratory rate, RRmand	1/min			01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Zeus IE	
7CH	1	Mandatory minute volume, MVmand	L/min			01.00	Perseus A500, Atlan, Zeus IE	

Code	Code page	Data description	Unit	LL	HL	New with version	Used by device	Remarks
7DH	1	Peak inspiratory pressure, PIP	mbar	Zeus IE, Fabius Family	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Zeus IE, Oxylog 3000 plus, Carina, Fabius Family	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300, Carina, Fabius Family	
7EH	1	Mandatory tidal volume, VTmand	L			01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Zeus IE	
7FH	1	Spontaneous tidal volume, VTspon	L			01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Zeus IE	
81H	1	Volume trapped in the lung by intrinsic PEEP, Vtrap	mL			01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800	

Code	Code page	Data description	Unit	LL	HL	New with version	Used by device	Remarks
83H	1	Mandatory expiratory tidal volume, VTemand	mL			01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800	
84H	1	Spontaneous expiratory tidal volume, VTespon	mL			01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800	
85H	1	Mandatory inspiratory tidal volume, VTimand	mL			01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800	
86H	1	Tidal volume during HFO, VThf	mL			01.00	Babylog VN500, Babylog VN600, Babylog VN800	

Code	Code page	Data description	Unit	LL	HL	New with version	Used by device	Remarks
88H	1	Tidal volume, VT	mL	Evita V300, Evita Infinity V500, Evita V600, Evita V800	Perseus A500, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Savina 300	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300, Carina, Fabius Family	
89H	1	Serial dead space volume, Vds	mL			01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800	
8AH	1	Ratio of serial dead space volume to expiratory tidal volume, Vds/VTe	%			01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800	
8BH	1	Inspiratory tidal volume, VTi (deprecated, will not be used by MEDIBUS.X devices)	mL			01.00	not used	
8CH	1	Non-invasive blood pressure - diastolic, NIBP D	mmHg	Zeus IE	Zeus IE	01.01	Zeus IE	
8DH	1	Negative inspiratory force, NIF	mbar			01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800	
91H	1	Non-invasive blood pressure - mean, NIBP M	mmHg	Zeus IE	Zeus IE	01.01	Zeus IE	
92H	1	Time since last non-invasive blood pressure measurement, Tnibp	min			01.01	Zeus IE	
96H	1	Non-invasive blood pressure - systolic, NIBP S	mmHg	Zeus IE	Zeus IE	01.01	Zeus IE	
9BH	1	Cardiac output, C.O.	L/min			01.01	Zeus IE	
9CH	1	Total halothane uptake, Hal tot (Uptake is calculated per anesthesia case.)	mL			01.01	Perseus A500, Atlan, Primus Family (added with SW 1.11 for Perseus A500)	

Code	Code page	Data description	Unit	LL	HL	New with version	Used by device	Remarks
9DH	1	Total desflurane uptake, Des tot (Uptake is calculated per anesthesia case.)	mL			01.01	Perseus A500, Atlan, Primus Family (added with SW 1.11 for Perseus A500)	
9EH	1	Total isoflurane uptake, Iso tot (Uptake is calculated per anesthesia case.)	mL			01.01	Perseus A500, Atlan, Primus Family (added with SW 1.11 for Perseus A500)	
9FH	1	Total sevoflurane uptake, Sev tot (Uptake is calculated per anesthesia case.)	mL			01.01	Perseus A500, Atlan, Primus Family (added with SW 1.11 for Perseus A500)	
A0H	1	Total enflurane uptake, Enf tot (Uptake is calculated per anesthesia case.)	mL			01.01	Perseus A500, Atlan, Primus Family (added with SW 1.11 for Perseus A500)	
A1H	1	Pulmonary arterial pressure - diastolic, PA D	mmHg	Zeus IE	Zeus IE	01.01	Zeus IE	
A2H	1	Duration of phototherapy, Dur PT The value is transmitted using 4 ASCII encoded digits HHMM. The two left most digits, i.e. HH, resemble the full hours of duration, the two right most digits, i.e. MM, resemble the duration in minutes since the last full hour duration. If duration is less than 1 hour, then the first two digits are both Ascii encoded '0'. Both fields (HH and MM) are padded with a leading, ASCII encoded '0', if necessary, i.e., if the duration is less than 10 hours, HH is zero padded and if the duration since the last full hour is less then 10 minutes, then MM is zero padded. A value of 9959 means duration took 99 hours and 59 minutes or longer.	h:min			01.10	BiliLux	
A3H	1	Irradiance during phototherapy, Irradiance PT	$\mu\text{W}/\text{cm}^2/\text{nm}$			01.10	BiliLux	
A5H	1	Central venous pressure, CVP	mmHg	Zeus IE	Zeus IE	01.01	Zeus IE	
A6H	1	Pulmonary arterial pressure - mean, PA M	mmHg	Zeus IE	Zeus IE	01.01	Zeus IE	
A7H	1	Right atrial pressure, RA	mmHg	Zeus IE	Zeus IE	01.01	Zeus IE	
ABH	1	Pulmonary arterial pressure - systolic, PA S	mmHg	Zeus IE	Zeus IE	01.01	Zeus IE	
ACH	1	MAC multiple derived from inspiratory concentrations, inxMAC	no unit			01.00	Perseus A500, Atlan, Primus Family	

Code	Code page	Data description	Unit	LL	HL	New with version	Used by device	Remarks
ADH	1	MAC multiple derived from expiratory concentrations, xMAC	no unit			01.00	Perseus A500, Atlan, Primus Family, Zeus IE	
AEH	1	Inspiratory desflurane concentration, inDes	%	Perseus A500, Atlan, Primus Family, Zeus IE	Perseus A500, Atlan, Primus Family, Zeus IE	01.00	Perseus A500, Atlan, Primus Family, Zeus IE, Fabius Family	
AFH	1	End-tidal desflurane concentration, etDes	%			01.00	Perseus A500, Atlan, Primus Family, Zeus IE, Fabius Family	
B0H	1	Inspiratory sevoflurane concentration, inSev	%	Perseus A500, Atlan, Primus Family, Zeus IE	Perseus A500, Atlan, Primus Family, Zeus IE	01.00	Perseus A500, Atlan, Primus Family, Zeus IE, Fabius Family	
B1H	1	End-tidal sevoflurane concentration, etSev	%			01.00	Perseus A500, Atlan, Primus Family, Zeus IE, Fabius Family	
B2H	1	Leakage of the breathing system	mL/min			01.00	Perseus A500, Atlan, Primus Family, Zeus IE	
B3H	1	Leakage minute volume in % of inspiratory minute volume, % leak	%			01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300	
B4H	1	Respiratory rate based on pressure measurement, RRP	1/min			01.00	Perseus A500, Atlan, Primus Family	

Code	Code page	Data description	Unit	LL	HL	New with version	Used by device	Remarks
B5H	1	Spontaneous respiratory rate, RR _{spon}	1/min		Oxylog 3000 plus	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Zeus IE, Oxylog 3000 plus	
B6H	1	Spontaneous breathing portion of minute volume, % MV _{spon}	%			01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800	
B7H	1	Spontaneous minute volume, MV _{spon}	L/min			01.00	Perseus A500, Atlan, Savina 300, Zeus IE, Oxylog 3000 plus (removed with SW 5.03 for Savina 300)	
B9H	1	Minute volume, MV	L/min	Perseus A500, Atlan, Primus Family, Zeus IE, Oxylog 3000 plus, Carina, Fabius Family	Perseus A500, Atlan, Primus Family, Zeus IE, Oxylog 3000 plus, Carina, Fabius Family	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300, Carina, Fabius Family	Edition 21: As of SW 5.03 Savina 300 does not send this code within responses to the commands 25H "Request current low Alarm Limits (Codepage 1)" and 26H Request current high Alarm Limits (Codepage 1)
BDH	1	Apnea duration	s			01.00	Perseus A500, Atlan, Primus Family	

Code	Code page	Data description	Unit	LL	HL	New with version	Used by device	Remarks
BEH	1	Temperature 2, T2	°C	Zeus IE	Zeus IE	01.01	Zeus IE	
C0H	1	Oxygen concentration, O2	%			01.08	Babyleo TN500, Isolette 8000 plus	
C1H	1	Airway temperature, T Airw	°C			01.00	Savina 300	
C2H	1	Skin temperature difference, ΔT Skin	°C			01.08	Babyleo TN500, Isolette 8000 plus, Babyroo TN300	
C3H	1	Temperature 1, T1	°C	Zeus IE	Zeus IE	01.01	Zeus IE	
C4H	1	Inspiratory/expiratory oxygen concentration difference, ΔO2	%			01.00	Perseus A500, Atlan, Primus Family	
C5H	1	Pulmonary wedge pressure, PWP	mmHg			01.01	Zeus IE	
C6H	1	Battery capacity	%			01.00	Savina 300, Carina	
C7H	1	Insp. O2 accuracy, O2bias	%			01.00	Savina 300	
C8H	1	Arterial diastolic pressure, ART D	mmHg	Zeus IE	Zeus IE	01.01	Zeus IE	
C9H	1	Rapid shallow breathing index, RSBI (synonymous parameter label: RSB)	1/min/L			01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300	Edition 19: New devices use a different terminology for the same parameter. Edition 6: The unit specification was missing parentheses. It has been corrected to highlight that RSB is a ratio of a frequency and a volume.
CDH	1	Arterial mean pressure, ART M	mmHg	Zeus IE	Zeus IE	01.01	Zeus IE	
CEH	1	Perfusion index, PI	%			01.13	Babyroo TN300	
D0H	1	Respiratory rate of triggered, mandatory breaths, RRtrig	1/min			01.06	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800 (added with SW 2.42 for Evita Infinity V500, Evita V300 and Babylog VN500)	
D2H	1	Arterial systolic pressure, ART S	mmHg	Zeus IE	Zeus IE	01.01	Zeus IE	

Code	Code page	Data description	Unit	LL	HL	New with version	Used by device	Remarks
D5H	1	Respiratory rate based on CO2 measurement, RRc	1/min			01.00	Perseus A500, Atlan, Primus Family, Zeus IE, Fabius Family	
D6H	1	Respiratory rate, RR	1/min		Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Carina	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300, Carina, Fabius Family (added with SW 1.10 for Zeus IE)	
D7H	1	Respiratory rate based on volume/flow measurement, RRf	1/min			01.00	Perseus A500, Atlan, Primus Family, Zeus IE	
D9H	1	Respiratory rate, RR (deprecated, will not be used by MEDIBUS.X devices)	1/min			01.00	not used	
DAH	1	Inspiratory CO2 concentration, inCO2	%		Perseus A500, Atlan, Primus Family, Zeus IE	01.00	Perseus A500, Atlan, Primus Family, Zeus IE, Fabius Family	

Code	Code page	Data description	Unit	LL	HL	New with version	Used by device	Remarks
DBH	1	End-tidal CO2 concentration, etCO2	%	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300, Fabius Family	
DCH	1	Heart rate, HR	1/min	Zeus IE	Zeus IE	01.01	Zeus IE	
DDH	1	N2O fresh-gas flow, FG V'N2O (Note: 9999 mL/min equals to 9999 mL/min or more.)	mL/min			01.00	Perseus A500, Atlan, Primus Family, Zeus IE, Fabius Family except Fabius plus, Fabius Plus XL	
DEH	1	Air fresh-gas flow, FG V'Air (Note: 9999 mL/min equals to 9999 mL/min or more.)	mL/min			01.00	Perseus A500, Atlan, Primus Family, Zeus IE, Fabius Family except Fabius plus, Fabius Plus XL	
E1H	1	Pulse Rate, PLS (value from channel 1, if multiple channels available)	1/min	Primus Family	Primus Family	01.01	Primus Family, Babyroo TN300	
E2H	1	Oxygen fresh-gas flow, FG V'O2 (Note: 9999 mL/min equals to 9999 mL/min or more.)	mL/min			01.00	Perseus A500, Atlan, Primus Family, Zeus IE, Fabius Family except Fabius plus, Fabius Plus XL	

Code	Code page	Data description	Unit	LL	HL	New with version	Used by device	Remarks
E3H	1	End-tidal CO2 concentration, etCO2	kPa	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300, Fabius Family	
E4H	1	Light level, Light	lux			01.09	Babyleo TN500	
E5H	1	Inspiratory CO2 concentration, inCO2	mmHg		Perseus A500, Atlan, Primus Family, Zeus IE	01.00	Perseus A500, Atlan, Primus Family, Zeus IE, Fabius Family	
E6H	1	End-tidal CO2 concentration, etCO2	mmHg	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300, Fabius Family	
E7H	1	Measured I : E I part	no unit			01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300	

Code	Code page	Data description	Unit	LL	HL	New with version	Used by device	Remarks
E8H	1	Measured I : E E part	no unit			01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300	
E9H	1	Inspiratory concentration of primary agent	%			01.00	Perseus A500, Atlan, Primus Family, Zeus IE	
EAH	1	End-tidal concentration of primary agent	%			01.00	Perseus A500, Atlan, Primus Family, Zeus IE	
EBH	1	Oxygen saturation from pulse oximetry, SpO2 (value from channel 1, if multiple channels available)	%	Primus Family, Zeus IE	Primus Family, Zeus IE	01.01	Primus Family, Zeus IE, Babyroo TN300	
ECH	1	Noise level, Noise	dB(A)			01.09	Babyleo TN500	
EDH	1	Inspiratory concentration of secondary agent	%			01.00	Perseus A500, Atlan, Primus Family, Zeus IE	
EEH	1	End-tidal concentration of secondary agent	%			01.00	Perseus A500, Atlan, Primus Family, Zeus IE	
EFH	1	End-tidal oxygen concentration, etO2	%			01.00	Perseus A500, Atlan, Primus Family, Zeus IE	

Code	Code page	Data description	Unit	LL	HL	New with version	Used by device	Remarks
F0H	1	Inspiratory oxygen fraction, FiO2	%	Perseus A500, Atlan, Primus Family, Zeus IE, Fabius Family	Perseus A500, Atlan, Primus Family, Zeus IE, Fabius Family	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Fabius Family	
F1H	1	Convective heater power, Heater conv	%			01.08	Babyleo TN500, Isolette 8000 plus	
F2H	1	Body weight, BW	g			01.08	Babyleo TN500, Isolette 8000 plus, Babyroo TN300	
F3H	1	Weight date (date, when latest weight measurement has been recorded) Four ASCII digit value: first two characters identify the day of month with leading '0' (30H), if necessary. Last two characters identify the month of the year with leading '0' (30H), if necessary.	DDMM			01.09	Babyleo TN500, Babyroo TN300	
F4H	1	Inspiratory halothane concentration, inHal	%	Perseus A500, Atlan, Primus Family, Zeus IE	Perseus A500, Atlan, Primus Family, Zeus IE	01.00	Perseus A500, Atlan, Primus Family, Zeus IE, Fabius Family	
F5H	1	End-tidal halothane concentration, etHal	%			01.00	Perseus A500, Atlan, Primus Family, Zeus IE, Fabius Family	
F6H	1	Inspiratory enflurane concentration, inEnf	%	Perseus A500, Atlan, Primus Family, Zeus IE	Perseus A500, Atlan, Primus Family, Zeus IE	01.00	Perseus A500, Atlan, Primus Family, Zeus IE, Fabius Family	

Code	Code page	Data description	Unit	LL	HL	New with version	Used by device	Remarks
F7H	1	End-tidal enflurane concentration, etEnf	%			01.00	Perseus A500, Atlan, Primus Family, Zeus IE, Fabius Family	
F8H	1	Inspiratory isoflurane concentration, inIso	%	Perseus A500, Atlan, Primus Family, Zeus IE	Perseus A500, Atlan, Primus Family, Zeus IE	01.00	Perseus A500, Atlan, Primus Family, Zeus IE, Fabius Family	
F9H	1	End-tidal isoflurane concentration, etIso	%			01.00	Perseus A500, Atlan, Primus Family, Zeus IE, Fabius Family	
FAH	1	Weight sample age	h			01.08	Isolette 8000 plus	
FBH	1	Inspiratory N2O concentration, inN2O	%	Zeus IE	Zeus IE	01.00	Perseus A500, Atlan, Primus Family, Zeus IE, Fabius Family	
FCH	1	End-tidal N2O concentration, etN2O	%			01.00	Perseus A500, Atlan, Primus Family, Zeus IE, Fabius Family	
FFH	1	Inspiratory CO2 concentration, inCO2	kPa		Perseus A500, Atlan, Primus Family, Zeus IE	01.00	Perseus A500, Atlan, Primus Family, Zeus IE, Fabius Family	
00H	2	Spontaneous tidal volume, VTspon	mL			01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800	
03H	2	Pressure amplitude during HFO, ΔPhf	mbar			01.00	Babylog VN500, Babylog VN600, Babylog VN800	

Code	Code page	Data description	Unit	LL	HL	New with version	Used by device	Remarks
04H	2	Inspiratory part of I:E ratio during spontaneous breathing, I (I:Espon)	no unit			01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800	Edition 19: Clarification of parameter semantics.
05H	2	Expiratory part of I:E ratio during spontaneous breathing, E (I:Espon)	no unit			01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800	Edition 19: Clarification of parameter semantics.
06H	2	Elastance, E	mbar/L			01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800	
07H	2	Time constant, TC	s			01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800	
09H	2	Expiratory time constant, TCe	s			01.01	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800	
0AH	2	Static compliance, Cstat	mL/mbar			01.00	Evita Infinity V500, Evita V600, Evita V800	

Code	Code page	Data description	Unit	LL	HL	New with version	Used by device	Remarks
0BH	2	Ratio of compliance during the last 20 % of inspiration over dynamic compliance, C20/Cdyn	no unit			01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800	
0CH	2	CO2 elimination coefficient during HFO, DCO2	10*mL^2/s			01.00	Babylog VN500, Babylog VN600, Babylog VN800	
0DH	2	Rapid shallow breathing index, RSBI (synonymous parameter label: RSB)	1/min/mL			01.02	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800 (added with SW 2.30 for Babylog VN500, added with SW 2.30 for Evita Infinity V500)	Edition 19: New devices use a different terminology for the same parameter. Edition 6: The unit specification was missing parentheses. It has been corrected to highlight that RSB is a ratio of a frequency and a volume.
0EH	2	Train-of-four ratio, TOF Ratio	%			01.01	Zeus IE	
0FH	2	Post-tetanic count, PTC	no unit			01.01	Zeus IE	
10H	2	NMT temperature, T NMT	°C			01.01	Zeus IE	
11H	2	Single twitch mode, Single	%			01.01	Zeus IE	
12H	2	Train-of-four count, TOF Cnt	no unit	Zeus IE	Zeus IE	01.01	Zeus IE	
13H	2	Bispectral index, BIS	no unit	Zeus IE	Zeus IE	01.01	Zeus IE	
14H	2	Electromyogram power, EMG	dB			01.01	Zeus IE	
15H	2	Signal quality index, SQI	%			01.01	Zeus IE	
16H	2	Burst suppression ratio (BIS), BSR	%			01.01	Zeus IE	
17H	2	Spectral edge frequency 95, SEF95	Hz			01.01	Zeus IE	
18H	2	Burst count (BIS), BCT	no unit			01.01	Zeus IE	
19H	2	Total power (BIS), EEG Pwr	dB			01.01	Zeus IE	
1AH	2	Goal for Δpressure support in SmartCare, SC-ΔPsupp goal	mbar			01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800	
1BH	2	Rated Δpressure support by internal controller in SmartCare, SC-ΔPsupp rated	mbar			01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800	

Code	Code page	Data description	Unit	LL	HL	New with version	Used by device	Remarks
1CH	2	Duration of session (hours) in SmartCare, SC-duration	h			01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800	
1DH	2	Duration of session (minutes) in SmartCare, SC-duration	min			01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800	
1EH	2	Spontaneous respiratory rate in SmartCare, SC-RRspon	1/min			01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800	
1FH	2	Tidal volume in SmartCare, SC-VT	mL			01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800	
21H	2	Expiratory tidal volume, VTe	mL	Evita V600, Evita V800, Babylog VN600, Babylog VN800		01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300	Edition 22: As of SW 3.00 Evita V600, Evita V800, Babylog VN600 and Babylog VN800 support this code also for the command 2CH "Request current low Alarm Limits (Codepage 2).
22H	2	Inspiratory tidal volume, VTi	mL	Evita V300, Evita Infinity V500, Evita V600, Evita V800	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Zeus IE	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Primus Family, Zeus IE, Carina	
23H	2	End-inspiratory pressure, EIP	mbar			01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800	

Code	Code page	Data description	Unit	LL	HL	New with version	Used by device	Remarks
24H	2	Effective expiratory time during APRV/AutoRelease, Tlow	s			01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800	
25H	2	CO2 production volume per breath, VTCO2	mL			01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800	
28H	2	End-tidal CO2 concentration in SmartCare, SC-etCO2	mmHg			01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800	
29H	2	ST deviation in lead I, STI	mm*0,1			01.01	Zeus IE	
2AH	2	ST deviation in lead II, STII	mm*0,1			01.01	Zeus IE	
2BH	2	ST deviation in lead III, STIII	mm*0,1			01.01	Zeus IE	
2CH	2	ST deviation in lead aVL, STaVL	mm*0,1			01.01	Zeus IE	
2DH	2	ST deviation in lead aVF, STaVF	mm*0,1			01.01	Zeus IE	
2EH	2	ST deviation in lead aVR, STaVR	mm*0,1			01.01	Zeus IE	
2FH	2	ST deviation in lead V, STV	mm*0,1			01.01	Zeus IE	
30H	2	ST deviation in lead V+, STV+	mm*0,1			01.01	Zeus IE	
31H	2	Upper pressure level during APRV, Phigh	mbar			01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800	
32H	2	Lower pressure level during APRV, Plow	mbar			01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800	
33H	2	Spontaneous mean tidal volume, VTspon mean	mL			01.00	Evita Infinity V500, Evita V600, Evita V800	
34H	2	Inspiratory spontaneous mean tidal volume, VTispon mean	mL			01.00	Evita Infinity V500, Evita V600, Evita V800	

Code	Code page	Data description	Unit	LL	HL	New with version	Used by device	Remarks
35H	2	Expiratory spontaneous mean Tidal Volume, VTespon mean	mL			01.00	Evita Infinity V500, Evita V600, Evita V800	
36H	2	Tidal volume per kg ideal body weight, VT/IBW (synonymous parameter label VT/Wt)	mL/kg			01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800	Edition 19: Clarification of parameter semantics. New devices use a different terminology for the same parameter.
37H	2	Leakage minute volume, MVleak	L/min			01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Carina (added with SW 5.03 for Savina 300)	
38H	2	Inspiratory time, Ti	s			01.00	Savina 300	
39H	2	Leakage of breathing system, bag and hoses in Man/Spont	mL/min			01.00	Perseus A500, Atlan, Primus Family	
3AH	2	Compliance of the breathing system including patient circuit, Cbs+hose	mL/hPa			01.00	Perseus A500, Atlan, Primus Family	Edition 6: Parameter name added to the description.
3BH	2	Compliance of the breathing hoses, Chose	mL/hPa			01.00	Perseus A500, Atlan	Edition 6: Parameter name added to the description.
3CH	2	Adaption Time, Tadapt	min			01.02	Carina	
3DH	2	Applied upper limit for the VT target range in Smart Ventilation Control, SVC-VT high	mL/kg			01.08	Zeus IE (added with SW 2.00 for Zeus IE)	
3EH	2	Tidal volume calculated by Smart Ventilation Control, SVC-VT	mL/kg			01.08	Zeus IE (added with SW 2.00 for Zeus IE)	
3FH	2	Applied lower limit for the VT target range in Smart Ventilation Control, SVC-VT low	mL/kg			01.08	Zeus IE (added with SW 2.00 for Zeus IE)	
40H	2	Applied upper limit for the etCO2 target range in Smart Ventilation Control, SVC-etCO2 high	mmHg			01.08	Zeus IE (added with SW 2.00 for Zeus IE)	

Code	Code page	Data description	Unit	LL	HL	New with version	Used by device	Remarks
41H	2	End-tidal carbon dioxide concentration calculated by Smart Ventilation Control, SVC-etCO2	mmHg			01.08	Zeus IE (added with SW 2.00 for Zeus IE)	
42H	2	Applied lower limit for the etCO2 target range in Smart Ventilation Control, SVC-etCO2 low	mmHg			01.08	Zeus IE (added with SW 2.00 for Zeus IE)	
43H	2	Total respiratory rate in Smart Ventilation Control, SVC-RRTotal	1/min			01.08	Zeus IE (added with SW 2.00 for Zeus IE)	
44H	2	Spontaneous respiratory rate set by Smart Ventilation Control, SVC-RRspon	1/min			01.08	Zeus IE (added with SW 2.00 for Zeus IE)	
45H	2	Relative pressure support above PEEP set by Smart Ventilation Control, SVC-ΔPsupp	mbar			01.08	Zeus IE (added with SW 2.00 for Zeus IE)	
46H	2	Inspiratory pressure set by Smart Ventilation Control, SVC-Plnsp	mbar			01.08	Zeus IE (added with SW 2.00 for Zeus IE)	
47H	2	Respiratory rate set by Smart Ventilation Control, SVC-RR	1/min			01.08	Zeus IE (added with SW 2.00 for Zeus IE)	
48H	2	Minimum respiratory rate set by Smart Ventilation Control, SVC-RRmin	1/min			01.08	Zeus IE (added with SW 2.00 for Zeus IE)	
49H	2	Inspiratory time set by Smart Ventilation Control, SVC-TI	s			01.08	Zeus IE (added with SW 2.00 for Zeus IE)	
4AH	2	Flow trigger set by Smart Ventilation Control, SVC-Trigger	L/min			01.08	Zeus IE (added with SW 2.00 for Zeus IE)	
4BH	2	currently active ventilation goal of Smart Ventilation Control, Ventilation goal Enumeration value with the following definition: 0: Controlled Ventilation, 1: Augmented Ventilation, 2: Encourage SB, 3: Prep. extubation	no Unit			01.08	Zeus IE (added with SW 2.00 for Zeus IE)	
4CH	2	Current state of Smart Ventilation Control, SVC state Enumeration value with the following definition: 0: Off, 1: On, 2: Suspended	no Unit			01.08	Zeus IE (added with SW 2.00 for Zeus IE)	

Code	Code page	Data description	Unit	LL	HL	New with version	Used by device	Remarks
4DH	2	VT classification by Smart Ventilation Control, SVC VT class. Enumeration value with the following definition: 0: Unknown, 1: very low, 2: low, 3: normal, 4: high, 5: very high	no Unit			01.08	Zeus IE (added with SW 2.00 for Zeus IE)	
4EH	2	EtCO2 classification by Smart Ventilation Control, SVC etCO2 class. Enumeration value with the following definition: 0: Unknown, 1: severe hypoventilated, 2: mild hypoventilated, 3: normoventilated, 4: mild hyperventilated, 5: severe hyperventilated	no Unit			01.08	Zeus IE (added with SW 2.00 for Zeus IE)	
4FH	2	Therapy phase of Smart Ventilation Control, SVC phase Enumeration value with the following definition: 0: Adaptation, 1: Ventilation Management, 2: Recovery	no Unit			01.08	Zeus IE (added with SW 2.00 for Zeus IE)	
50H	2	Ventilation mode used by Smart Ventilation Control, SVC mode Enumeration value with the following definition: 0: PC-BIPAP/PS, 1: SPN-CPAP/ PS	no Unit			01.08	Zeus IE (added with SW 2.00 for Zeus IE)	
51H	2	Spontaneous breathing detection of Smart Ventilation Control, SVC spont. breathing Enumeration value with the following definition: 0: No, 1: Yes	no Unit			01.08	Zeus IE (added with SW 2.00 for Zeus IE)	
52H	2	Resistance of the breathing hoses, Rhose	mbar/(L/s)			01.03	Perseus A500 (added with SW 1.11 for Perseus A500)	Edition 15: Specification of unit now contains brackets.
53H	2	Difference between inspired and expired tidal volume (VTi - VT _e), ΔVT	mL			01.03	Perseus A500, Atlan (added with SW 1.11 for Perseus A500)	

Real-time data

Code	Data description	Unit	New with version	Used by device	Remarks
00H	Airway pressure, Paw	mbar	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300, Carina, Fabius Family	
01H	Flow (inspiratory/expiratory), Flow	L/min	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300, Carina, Fabius Family (added with SW 1.10 for Oxylog VE300)	
02H	Oxygen saturation pulse, Pleth (value from channel 1, if multiple channels available)	no unit	01.01	Primus Family, Babyroo TN300	

Code	Data description	Unit	New with version	Used by device	Remarks
00H	Airway pressure, Paw	mbar	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300, Carina, Fabiuss Family	
03H	Respiratory volume since start of inspiration, Volume	mL	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300	
05H	Oxygen concentration (inspiratory/expiratory), O2	%	01.00	Perseus A500, Atlan, Primus Family, Zeus IE	
06H	CO2 concentration, CO2	mmHg	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300	

Code	Data description	Unit	New with version	Used by device	Remarks
00H	Airway pressure, Paw	mbar	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300, Carina, Fabius Family	
07H	CO2 concentration, CO2	kPa	01.00	Perseus A500, Atlan, Evita V600, Evita V800, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300 (added with SW 2.00 for Evita V800, Evita V600, Babylog VN800 and Babylog VN600)	

Code	Data description	Unit	New with version	Used by device	Remarks
00H	Airway pressure, Paw	mbar	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300, Carina, Fabius Family	
08H	CO2 concentration, CO2	%	01.00	Perseus A500, Atlan, Evita V600, Evita V800, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300 (added with SW 2.00 for Evita V800, Evita V600, Babylog VN800 and Babylog VN600)	
0AH	Concentration of primary agent (inspiratory/expiratory)	%	01.00	Perseus A500, Atlan, Primus Family	
0BH	Halothane concentration (inspiratory/expiratory), Hal	%	01.00	Perseus A500, Atlan, Primus Family	
0CH	Enflurane concentration (inspiratory/expiratory), Enf	%	01.00	Perseus A500, Atlan, Primus Family	
0DH	Isoflurane concentration (inspiratory/expiratory), Iso	%	01.00	Perseus A500, Atlan, Primus Family	

Code	Data description	Unit	New with version	Used by device	Remarks
00H	Airway pressure, Paw	mbar	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300, Carina, Fabius Family	
0EH	Desflurane concentration (inspiratory/expiratory), Des	%	01.00	Perseus A500, Atlan, Primus Family	
0FH	Sevoflurane concentration (inspiratory/expiratory), Sev	%	01.00	Perseus A500, Atlan, Primus Family	
1CH	Tracheal pressure, Ptrach	mbar	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300 (added with SW 5.00 for Savina 300)	
1EH	Inspiratory device flow, Flow dev insp	L/min	01.04	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800 (added with SW 2.40 for Babylog VN500, added with SW 2.40 for Evita Infinity V500 and Evita V300)	

Code	Data description	Unit	New with version	Used by device	Remarks
00H	Airway pressure, Paw	mbar	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300, Carina, Fabius Family	
2AH	Concentration of primary agent (inspiratory/expiratory)	kPa	01.00	Perseus A500, Atlan, Primus Family	
2BH	Halothane concentration (inspiratory/expiratory), Hal	kPa	01.00	Perseus A500, Atlan, Primus Family	
2CH	Enflurane concentration (inspiratory/expiratory), Enf	kPa	01.00	Perseus A500, Atlan, Primus Family	
2DH	Isoflurane concentration (inspiratory/expiratory), Iso	kPa	01.00	Perseus A500, Atlan, Primus Family	
2EH	Desflurane concentration (inspiratory/expiratory), Des	kPa	01.00	Perseus A500, Atlan, Primus Family	
2FH	Sevoflurane concentration (inspiratory/expiratory), Sev	kPa	01.00	Perseus A500, Atlan, Primus Family	
C1H	Enable/Disable datastream 1 to 4		01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800	

Code	Data description	Unit	New with version	Used by device	Remarks
00H	Airway pressure, Paw	mbar	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300, Carina, Fabius Family	
C2H	Enable/Disable datastream 5 to 8		01.03	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800	Edition 6: Even though this sync command was not documented with the MEDIBUS.X profile before version 1.03, Evita Infinity V500 and Babylog VN500 already used it in their implementation of the MEDIBUS realtime-extension
C6 / C0H	Sync signal start of ventilator inspiratory cycle (command code / argument)		01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300, Carina, Fabius Family (added with SW 2.00 for Zeus IE)	
C6 / C1H	Sync signal start of ventilator expiratory cycle (command code / argument)		01.01	Oxylog 3000 plus, Oxylog VE300, Carina	Edition 6: Evita Infinity V500 and Babylog VN500 do not send this sync signal

Alarm status response

Code	Code page	Alarm description	Alarm phrase	New with version	Used by device	Remarks
00H	1	Apnea by different or undefined sources	APNEA	01.00	Perseus A500, Atlan, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300	
01H	1	No SpO2 pulse	NO SPO2 PULS	01.01	Primus Family, Zeus IE	
02H	1	SpO2 pulse < low limit	SPO2 PULS LO	01.01	Primus Family, Babyroo TN300	
03H	1	General systolic pressure (channel 1) > high limit	SYS GP1 HIGH	01.01	Zeus IE	
04H	1	General systolic pressure (channel 1) < low limit	SYS GP1 LOW	01.01	Zeus IE	
05H	1	Heart rate / pulse < low limit	HRT RATE LOW	01.01	Zeus IE	
06H	1	Pulmonary arterial pressure - systolic > high limit	SYS PA HI	01.01	Zeus IE	
07H	1	SpO2 < low limit	SPO2 LOW	01.01	Primus Family, Zeus IE, Babyroo TN300	
08H	1	Inspiratory oxygen concentration < low limit	% O2 LOW	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Fabius Family	
09H	1	Inspiratory halothane concentration > high limit	% HAL HIGH	01.00	Perseus A500, Atlan, Primus Family, Zeus IE	
0AH	1	Pulmonary arterial pressure - diastolic > high limit	DIAS PA HI	01.01	Zeus IE	

Code	Code page	Alarm description	Alarm phrase	New with version	Used by device	Remarks
0BH	1	Inspiratory enflurane concentration > high limit	% ENF HIGH	01.00	Perseus A500, Atlan, Primus Family, Zeus IE	
0CH	1	Inspiratory isoflurane > high limit	% ISO HIGH	01.00	Perseus A500, Atlan, Primus Family, Zeus IE	
0DH	1	Apnea - detected by CO2 monitoring.	APNEA CO2	01.00	Perseus A500, Atlan, Primus Family, Zeus IE	
0EH	1	Apnea - detected by volume monitoring	APNEA VOL	01.00	Perseus A500, Atlan, Primus Family, Zeus IE, Fabius Family	
0FH	1	Apnea - detected by pressure monitoring.	APNEA PRES	01.00	Perseus A500, Atlan, Primus Family, Zeus IE, Fabius Family	
10H	1	Airway pressure > high limit	PAW HIGH	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300, Carina, Fabius Family	
11H	1	Check fresh-gas setting	FRESH GAS ?	01.00	Perseus A500, Atlan, Primus Family, Zeus IE, Fabius Family	

Code	Code page	Alarm description	Alarm phrase	New with version	Used by device	Remarks
12H	1	Check air supply	AIR SUPPLY ?	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Primus Family, Zeus IE	
13H	1	Check O2 supply	O2 SUPPLY ?	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300, Fabius Family	
14H	1	Pulmonary arterial pressure - mean > high limit	MEAN PA HI	01.01	Zeus IE	
15H	1	Pulmonary arterial pressure - systolic < low limit	SYS PA LOW	01.01	Zeus IE	
16H	1	Pulmonary arterial pressure - diastolic < low limit	DIAS PA LO	01.01	Zeus IE	
17H	1	Pulmonary arterial pressure - mean < low limit	MEAN PA LO	01.01	Zeus IE	
18H	1	Expiratory pressure > high limit	PRESS EXP HI	01.05	Fabius Family	

Code	Code page	Alarm description	Alarm phrase	New with version	Used by device	Remarks
19H	1	Minute volume < low limit	MIN VOL LOW	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300, Carina, Fabius Family	
1AH	1	NIBP systolic < low limit	NBP SYS LOW	01.01	Zeus IE	
1BH	1	NIBP diastolic < low limit	NBP DIA LOW	01.01	Zeus IE	
1CH	1	Arterial systolic pressure < low limit	SYS ART LOW	01.01	Zeus IE	
1DH	1	General diastolic pressure (channel 1) > high limit	DIAS GP1 HI	01.01	Zeus IE	
1EH	1	SpO2 pulse > high limit	SPO2 PULS HI	01.01	Primus Family, Babyroo TN300	
1FH	1	Inspiratory sevoflurane concentration > high limit	% SEV HIGH	01.00	Perseus A500, Atlan, Primus Family, Zeus IE	
20H	1	General diastolic pressure (channel 1) < low limit	DIAS GP1 LOW	01.01	Zeus IE	
21H	1	Heart rate > high limit	HRT RATE HI	01.01	Zeus IE	
22H	1	SpO2 > high limit	SPO2 HIGH	01.01	Primus Family, Zeus IE, Babyroo TN300	
23H	1	NIBP systolic > high limit	NBP SYS HI	01.01	Zeus IE	
24H	1	Inspiratory desflurane concentration > high limit	% DES HIGH	01.00	Perseus A500, Atlan, Primus Family, Zeus IE	
25H	1	Arterial systolic pressure > high limit	SYS ART HI	01.01	Zeus IE	
26H	1	General mean pressure (channel 1) > high limit	MEAN GP1 HI	01.01	Zeus IE	

Code	Code page	Alarm description	Alarm phrase	New with version	Used by device	Remarks
27H	1	End-tidal CO2 < low limit	ET CO2 LOW	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300	
28H	1	End-tidal CO2 > high limit	ET CO2 HIGH	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300	
29H	1	Inspiratory halothane concentration < low limit	% HAL LOW	01.00	Perseus A500, Atlan, Primus Family, Zeus IE	
2AH	1	Inspiratory enflurane concentration < low limit	% ENF LOW	01.00	Perseus A500, Atlan, Primus Family, Zeus IE	
2BH	1	Inspiratory isoflurane < low limit	% ISO LOW	01.00	Perseus A500, Atlan, Primus Family, Zeus IE	
2CH	1	Inspiratory desflurane < low limit	% DES LOW	01.00	Perseus A500, Atlan, Primus Family, Zeus IE	
2DH	1	NIBP diastolic > high limit	NBP DIA HI	01.01	Zeus IE	

Code	Code page	Alarm description	Alarm phrase	New with version	Used by device	Remarks
2EH	1	Temperature 1 < low limit	TEMP1 LOW	01.01	Zeus IE	
2FH	1	Temperature 2 < low limit	TEMP2 LOW	01.01	Zeus IE	
30H	1	Air temperature high	AIR TEMP HI	01.08	Babyleo TN500, Isolette 8000 plus	
31H	1	General mean pressure (channel 1) < low limit	MEAN GP1 LOW	01.01	Zeus IE	
32H	1	Inspiratory sevoflurane < low limit	% SEV LOW	01.00	Perseus A500, Atlan, Primus Family, Zeus IE	
33H	1	Volume not constant	VOL INCONST	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Oxylog VE300	
35H	1	SpO2 sensor disconnected, not supported or fault	SPO2SEN DISC	01.01	Primus Family, Zeus IE, Babyroo TN300	
36H	1	CO2 Filter check failed	CO2 FIL. ERR	01.01	Oxylog 3000 plus, Oxylog VE300	
37H	1	Inspiratory oxygen concentration > high limit	% O2 HIGH	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Fabius Family	
38H	1	General systolic pressure (channel 2) > high limit	SYS GP2 HIGH	01.01	Zeus IE	
39H	1	General systolic pressure (channel 2) < low limit	SYS GP2 LOW	01.01	Zeus IE	
3AH	1	Pressure support terminated by time	PS LIMITED	01.00	Savina 300, Zeus IE	
3BH	1	Oxygen analyzer not calibrated	CAL % O2 ?	01.01	Carina, Fabius Family	

Code	Code page	Alarm description	Alarm phrase	New with version	Used by device	Remarks
3CH	1	Inspiratory CO2 > high limit	INSP CO2 HI	01.00	Perseus A500, Atlan, Primus Family, Zeus IE	
3DH	1	CO2 patient sensor line blocked	CO2 LINE BLK	01.00	Perseus A500, Atlan, Primus Family, Zeus IE	
3EH	1	CO2 zero-calibration requested	CO2 ZERO CAL	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Oxylog 3000 plus, Oxylog VE300	
3FH	1	NIBP mean < low limit	NBP MEAN LO	01.01	Zeus IE	
40H	1	General diastolic pressure (channel 2) > high limit	DIAS GP2 HI	01.01	Zeus IE	
41H	1	NIBP mean > high limit	NBP MEAN HI	01.01	Zeus IE	
42H	1	Check flow sensor	FLOW SENSOR?	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Savina 300, Zeus IE	
43H	1	O2 sensor inoperable	O2 SENS ERR	01.08	Babyleo TN500, Isolette 8000 plus	
44H	1	General diastolic pressure (channel 2) < low limit	DIAS GP2 LOW	01.01	Zeus IE	
45H	1	General mean pressure (channel 2) > high limit	MEAN GP2 HI	01.01	Zeus IE	
46H	1	Temperature 1 failure	TEMP1 ERR	01.01	Zeus IE	
47H	1	Temperature 2 failure	TEMP2 ERR	01.01	Zeus IE	
48H	1	Air temperature measurement failure	AIR TEMP ERR	01.08	Isolette 8000 plus	
49H	1	Check NIBP cuff size!	BP CUFF ERR	01.01	Zeus IE	
4AH	1	Unable to take NiBP measurement	BP CUFF DISC	01.01	Zeus IE	

Code	Code page	Alarm description	Alarm phrase	New with version	Used by device	Remarks
4BH	1	Battery capacity low	BATTERY LOW	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300, Carina, Fabius Family	
4CH	1	General mean pressure (channel 2) < low limit	MEAN GP2 LOW	01.01	Zeus IE	
4DH	1	Left ventricular pressure - systolic > high limit	SYS LV HIGH	01.01	Zeus IE	
4EH	1	Left ventricular pressure - systolic < low limit	SYS LV LOW	01.01	Zeus IE	
4FH	1	LL electrode failure	ECG ELEC. LL	01.01	Zeus IE	
50H	1	LA electrode failure	ECG ELEC. LA	01.01	Zeus IE	
51H	1	RA electrode failure	ECG ELEC. RA	01.01	Zeus IE	
52H	1	Left ventricular pressure - diastolic > high limit	DIAS LV HIGH	01.01	Zeus IE	
53H	1	ECG S-T analysis impossible	NO S-T	01.01	Zeus IE	
54H	1	Central venous pressure < low limit	CVP LO	01.01	Zeus IE	
55H	1	Central venous pressure > high limit	CVP HI	01.01	Zeus IE	
56H	1	Left ventricular pressure - diastolic < low limit	DIAS LV LOW	01.01	Zeus IE	
57H	1	CO2 alarm disabled	CO2 ALRM OFF	01.01	Primus Family, Zeus IE	
58H	1	No agent compensation	NO_AGAS_COMP	01.13	Atlan (added with SW 2.00 for Atlan)	
59H	1	No carrier-gas compensation	NO_CGAS_COMP	01.14	Atlan (added with SW 2.10 for Atlan)	
5AH	1	Cardiovasc. alarms disabled alarms off	HEMO ALRM OF	01.01	Zeus IE	
5BH	1	Pulse oximeter alarm disabled	SPO2 ALRM OF	01.01	Primus Family	
5CH	1	Left ventricular pressure - mean > high limit	MEAN LV HIGH	01.01	Zeus IE	
5EH	1	Volume alarm disabled	VOL ALRM OFF	01.05	Fabius Family	
60H	1	NIBP hose unplugged	NBP STANDBY	01.01	Zeus IE	

Code	Code page	Alarm description	Alarm phrase	New with version	Used by device	Remarks
61H	1	Left ventricular pressure - mean < low limit	MEAN LV LOW	01.01	Zeus IE	
62H	1	Left atrial pressure > high limit	LA HIGH	01.01	Zeus IE	
63H	1	Left atrial pressure < low limit	LA LOW	01.01	Zeus IE	
64H	1	Clean CO2 sensor (window occluded)	CLEAN CO2	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Oxylog 3000 plus, Oxylog VE300	
65H	1	Speaker failure	SPEAKER FAIL	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300, Carina, Babyleo TN500, Babyroo TN300, Fabius Family	
66H	1	Mixed agent detected	MIXED AGENT	01.00	Perseus A500, Atlan, Primus Family, Zeus IE	
67H	1	Multi-gas monitor device failure	GAS MON ERR	01.00	Perseus A500, Atlan, Primus Family	
68H	1	SpO2 device failure	SPO2 ERR	01.01	Primus Family, Zeus IE, Babyroo TN300	

Code	Code page	Alarm description	Alarm phrase	New with version	Used by device	Remarks
69H	1	N2O-measurement inoperable	N2O ERR	01.01	Perseus A500, Atlan, Primus Family, Zeus IE (added with SW 2.00 for Perseus A500)	
6AH	1	CO2 device failure	CO2 ERR	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300 (added with SW 2.00 for Perseus A500)	
6BH	1	Air temperature low	AIR TEMP LOW	01.08	Babyleo TN500, Isolette 8000 plus	
6CH	1	NiBP device failure	NBP ERR	01.01	Zeus IE	
6EH	1	Agent-measurement inoperable	AGT ERR	01.01	Perseus A500, Atlan, Primus Family, Zeus IE (added with SW 2.00 for Perseus A500)	
6FH	1	Right ventricular pressure - systolic > high limit	SYS RV HIGH	01.01	Zeus IE	
70H	1	Oxygen concentration low	OXYGEN LOW	01.08	Babyleo TN500, Isolette 8000 plus	
71H	1	Right ventricular pressure - systolic < low limit	SYS RV LOW	01.01	Zeus IE	
72H	1	Right ventricular pressure - diastolic > high limit	DIAS RV HIGH	01.01	Zeus IE	
73H	1	Right ventricular pressure - diastolic < low limit	DIAS RV LOW	01.01	Zeus IE	
74H	1	Right ventricular pressure - mean > high limit	MEAN RV HIGH	01.01	Zeus IE	
75H	1	Right ventricular pressure - mean < low limit	MEAN RV LOW	01.01	Zeus IE	
76H	1	Intracranial pressure > high limit	ICP HIGH	01.01	Zeus IE	
77H	1	Intracranial pressure < low limit	ICP LOW	01.01	Zeus IE	

Code	Code page	Alarm description	Alarm phrase	New with version	Used by device	Remarks
78H	1	Communication error RS232 port	RS232COM ERR	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Primus Family, Fabius Family	
79H	1	Cerebral perfusion pressure > high limit	CPP HIGH	01.01	Zeus IE	
7AH	1	Cerebral perfusion pressure < low limit	CPP LOW	01.01	Zeus IE	
7BH	1	Hemomed Pod 1 failure	HEMO POD ERR	01.01	Zeus IE	
7CH	1	Internal system communication error	INT COM ERR	01.01	Zeus IE	
7DH	1	Arterial pressure no pulse detected	NO ART PULS	01.01	Zeus IE	
7FH	1	Problems with intracranial pressure measurement	ICP ERR	01.01	Zeus IE	
80H	1	Problems with left atrial pressure measurement	LA ERR	01.01	Zeus IE	
81H	1	Problems with right ventricular pressure measurement	RV ERR	01.01	Zeus IE	
82H	1	Problems with left ventricular pressure measurement	LV ERR	01.01	Zeus IE	
83H	1	Problems with general blood pressure measurement (channel 1)	GP1 ERR	01.01	Zeus IE	
84H	1	Problems with general blood pressure measurement (channel 2)	GP2 ERR	01.01	Zeus IE	
85H	1	No pulse detected for right ventricular pressure measurement	NO RV PULS	01.01	Zeus IE	
86H	1	Arterial mean pressure < low limit	MEAN ART LO	01.01	Zeus IE	
87H	1	Arterial mean pressure > high limit	MEAN ART HI	01.01	Zeus IE	
88H	1	Arterial diastolic pressure < low limit	DIAS ART LO	01.01	Zeus IE	
89H	1	No pulse detected for left ventricular pressure measurement	NO LV PULS	01.01	Zeus IE	
8AH	1	Arterial diastolic pressure > high limit	DIAS ART HI	01.01	Zeus IE	
8BH	1	No pulse detected for general blood pressure measurement (channel 1)	NO GP1 PULS	01.01	Zeus IE	
8CH	1	No pulse detected for general blood pressure measurement (channel 2)	NO GP2 PULS	01.01	Zeus IE	

Code	Code page	Alarm description	Alarm phrase	New with version	Used by device	Remarks
90H	1	Respiratory rate > high limit	RESP RATE HI	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Oxylog 3000 plus, Oxylog VE300, Carina	
92H	1	ECG leads failure	ECG ELECTROD	01.01	Zeus IE	
94H	1	S-T segment high alarm	S-T HIGH	01.01	Zeus IE	
95H	1	Rescue ventilation active	RESCUE VENT	01.01	Carina	
96H	1	Asystoly	ASYSTOLE	01.01	Zeus IE	
98H	1	Apnea detected by respiratory monitoring	APNEA RESP	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Carina	
99H	1	S-T segment low alarm	S-T LOW	01.01	Zeus IE	
9AH	1	Disconnection ventilator / airway pressure < low limit	PAW LOW	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Zeus IE, Oxylog 3000 plus, Oxylog VE300, Carina (added with SW 2.00 for Atlan)	

Code	Code page	Alarm description	Alarm phrase	New with version	Used by device	Remarks
9BH	1	Minute volume > high limit	MIN VOL HIGH	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300, Carina, Fabius Family	
9EH	1	Water level low	WATER LOW	01.09	Babyleo TN500	
9FH	1	Problems with respirator	VENT ERR	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300, Carina, Fabius Family	
A0H	1	Ventilator Communication lost	COM VENT ERR	01.01	Primus Family, Zeus IE	
A1H	1	Inspiratory O2 measurement inoperable	% O2 ERR	01.00	Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300	

Code	Code page	Alarm description	Alarm phrase	New with version	Used by device	Remarks
A2H	1	Flow calibration necessary	VOL CAL ?	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Zeus IE, Fabius Family	
A3H	1	Airway pressure negative	PAW NEGATIVE	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Primus Family, Zeus IE, Fabius Family	
A4H	1	Volume measurement inoperable	VOL ERR	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300	
A6H	1	Check water supply	CHECK WATER	01.08	Isolette 8000 plus	
A7H	1	AutoAdapt active	AUTO ADAPT	01.02	Carina	
A8H	1	Oxygen concentration high	OXYGEN HIGH	01.08	Babyleo TN500, Isolette 8000 plus	
A9H	1	Leakage valve blocked	LEAK VALVE?	01.02	Carina	
ABH	1	Convective heater failure	AIR HTR INOP	01.08	Babyleo TN500, Isolette 8000 plus	
ACH	1	Minute volume alarm disabled	MV ALARM OFF	01.00	Savina 300, Primus Family, Zeus IE	

Code	Code page	Alarm description	Alarm phrase	New with version	Used by device	Remarks
ADH	1	Pressure measurement inoperable	PRESS ERR	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300, Carina, Fabius Family	
AEH	1	Continuous nebulization active	CONT. NEBUL.	01.04	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800 (added with SW 2.40 for Babylog VN500, added with SW 2.40 for Evita Infinity V500 and Evita V300)	
AFH	1	Check watertrap	WATER TRAP ?	01.00	Perseus A500, Atlan	
B0H	1	Check expiratory valve	EXP-VALVE ?	01.00	Perseus A500, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Savina 300, Zeus IE, Carina, Fabius Family	Edition 20: Babylog VN500, Babylog VN600 and Babylog VN800 do not send this code
B1H	1	Hood open - check hood	CHECK HOOD	01.09	Babyleo TN500	
B2H	1	Weaning temperature not reached	CHECK WEAN	01.09	Babyleo TN500	

Code	Code page	Alarm description	Alarm phrase	New with version	Used by device	Remarks
B3H	1	Gas calibration in progress	GAS MON CAL	01.11	Perseus A500, Atlan, Evita V600, Evita V800, Babylog VN600, Babylog VN800, Zeus IE, Oxylog VE300 (added with SW 1.01 for Atlan, added with SW 2.03 for Perseus A500, added with SW 2.02 for Zeus IE, added with SW 2.00 for Oxylog VE300)	
B4H	1	Heated mattress faulty	MAT HTR INOP	01.09	Babyleo TN500, Babyroo TN300	
B5H	1	Sensor module failure	SENS MOD ERR	01.08	Babyleo TN500, Isolette 8000 plus	
B6H	1	Warm-up procedure completed	WARM-UP DONE	01.09	Babyleo TN500, Babyroo TN300	
B7H	1	Temperature of ventilator blower too high	VENT TEMP HI	01.00	Perseus A500, Zeus IE	
B8H	1	Airway temperature measurement inop	AW-TEMP INOP	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300	
B9H	1	Check airway temperature sensor	AW-TEMP SENS	01.00	Evita V300, Evita Infinity V500, Babylog VN500	

Code	Code page	Alarm description	Alarm phrase	New with version	Used by device	Remarks
BAH	1	Airway temperature > high limit	AW-TEMP HIGH	01.00	Perseus A500, Atlan, Evita V600, Evita V800, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Carina (added with SW 1.11 for Perseus A500, added with SW 2.00 for Evita V800, Evita V600, Babylog VN800 and Babylog VN600)	
BBH	1	Session cannot be continued. Patient values do not match the guideline.	SCAVNOTADAPT	01.08	Zeus IE (added with SW 2.00 for Zeus IE)	
BCH	1	Breathing System not locked	SYSTEM OPEN	01.01	Zeus IE	
BDH	1	Check O2 supply (deprecated, will not be used by MEDIBUS.X devices)	LO O2 SUPPLY	01.00	not used	
BEH	1	O2 measurement inoperable	% O2 ERR	01.00	Perseus A500, Atlan, Primus Family, Zeus IE, Fabius Family	
BFH	1	Weaning procedure failed	WEANING FAIL	01.09	Babyleo TN500	
C0H	1	Inspiratory O2 < low limit	% O2 LOW	01.01	Zeus IE, Carina	
C1H	1	Flow measurement inoperable	VOL ERR	01.00	Perseus A500, Atlan, Primus Family, Zeus IE, Fabius Family	
C2H	1	Gas mixing and/or delivery failure	MIXER INOP	01.00	Perseus A500, Atlan, Zeus IE	

Code	Code page	Alarm description	Alarm phrase	New with version	Used by device	Remarks
C4H	1	Pressure limited respiratory volume	PRESSURE LIM	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Primus Family, Zeus IE, Fabius Family	
C5H	1	Weaning procedure completed	WEANING DONE	01.09	Babyleo TN500	
C7H	1	Fail to cycle	CYCLE FAILED	01.00	Savina 300	
C8H	1	Gas mixer inoperable	MIXER INOP	01.01	Primus Family, Oxylog 3000 plus	
C9H	1	Internal device temperature too high	INT.TMP.HIGH	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Primus Family, Carina, Babyleo TN500, Babyroo TN300	
CAH	1	Problems with fan	FAN ERR	01.00	Perseus A500, Atlan, Savina 300, Primus Family, Zeus IE, Babyleo TN500, Isolette 8000 plus	
CBH	1	Temperature 1 > high limit	TEMP1 HIGH	01.01	Zeus IE	
CCH	1	Temperature 2 > high limit	TEMP2 HIGH	01.01	Zeus IE	
CDH	1	O2 measurement temporarily unavailable	O2 UNAVAIL.	01.10	Perseus A500, Atlan (added with SW 2.00 for Perseus A500)	
CEH	1	N2O measurement temporarily unavailable	N2O UNAVAIL.	01.10	Perseus A500, Atlan (added with SW 2.00 for Perseus A500)	

Code	Code page	Alarm description	Alarm phrase	New with version	Used by device	Remarks
CFH	1	Agent measurement temporarily unavailable	AGT UNAVAIL.	01.10	Perseus A500, Atlan (added with SW 2.00 for Perseus A500)	
D0H	1	V-electrode failure	ECG ELEC. V	01.01	Zeus IE	
D1H	1	Patient gas measurement overflow	GASMON.+++	01.10	Perseus A500, Atlan (added with SW 2.00 for Perseus A500)	
D2H	1	VFIB / VTACH - arrhythmia detected	VENTR. FIBR.	01.01	Zeus IE	
D3H	1	Patient misconnected? Patient activity detected at inactive breathing system	PAT MISCONN?	01.11	Atlan	Edition 21: This code was introduced with version 01.11 of MEDIBUS.X. In former editions the documentation was incomplete with that regard.
D5H	1	Multiple gas measurement values can temporarily not be measured with accuracy complying to ISO standard.	GASMON.INACC	01.07	Perseus A500, Atlan, Primus Family, Zeus IE (added with SW 2.00 for Perseus A500, added with SW 2.00 for Zeus IE, added with SW 4.52 for the Primus Family)	
D6H	1	Problems with arterial pressure measurement	ART ERR	01.01	Zeus IE	
D7H	1	Problems with central venous pressure measurement	CVP ERR	01.01	Zeus IE	
D8H	1	Agent measurement values can temporarily not be measured with accuracy complying to ISO standard.	AGENT INACC.	01.07	Perseus A500, Atlan, Primus Family, Zeus IE (added with SW 2.00 for Perseus A500, added with SW 2.00 for Zeus IE, added with SW 4.52 for the Primus Family)	

Code	Code page	Alarm description	Alarm phrase	New with version	Used by device	Remarks
D9H	1	CO2 sensor disconnected or fault	CO2 SENSOR ?	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Oxylog 3000 plus, Oxylog VE300	
DAH	1	PEEP > high limit	PEEP HIGH	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Fabius Family (removed with SW 2.00 for Atlan)	
DBH	1	SpO2 < desaturation limit	SPO2 DESAT	01.13	Babyroo TN300	
DDH	1	Cardiac output device failure	CO ERR	01.01	Zeus IE	
DEH	1	Cardiac output check injectate probe	INJ. PROBE ?	01.01	Zeus IE	
DFH	1	Cardiac output transducer fail	CATHETER ?	01.01	Zeus IE	
E0H	1	SpO2 sensor has detected interference	SPO2 INTERFE	01.13	Babyroo TN300	
E1H	1	SpO2 accessory close to expiry	SPO2 EXPSOON	01.13	Babyroo TN300	
E2H	1	Problems with pulmonary arterial pressure measurement	PA ERR	01.01	Zeus IE	
E3H	1	Pulmonary arterial pressure no pulse detected	NO PA PULS	01.01	Zeus IE	

Code	Code page	Alarm description	Alarm phrase	New with version	Used by device	Remarks
E4H	1	N2O measurement values can temporarily not be measured with accuracy complying to ISO standard.	N2O INACC.	01.07	Perseus A500, Atlan, Primus Family, Zeus IE (added with SW 2.00 for Perseus A500, added with SW 2.00 for Zeus IE, added with SW 4.52 for the Primus Family)	
E5H	1	O2 measurement values can temporarily not be measured with accuracy complying to ISO standard.	O2 INACC.	01.07	Perseus A500, Atlan, Evita V600, Evita V800, Babylog VN600, Babylog VN800, Primus Family, Zeus IE (added with SW 2.00 for Perseus A500, added with SW 2.00 for Zeus IE, added with SW 4.52 for the Primus Family)	
E6H	1	Air supply pressure > high limit	AIR PRESS HI	01.00	Perseus A500, Atlan, Babylog VN500, Babylog VN600, Babylog VN800, Zeus IE	Edition 20: Evita V300, Evita Infinity V500, Evita V600 and Evita V800 do not send this code.
E7H	1	High O2 supply pressure	HI O2 SUPPLY	01.00	Perseus A500, Atlan, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Zeus IE, Carina	Edition 20: Evita V300, Evita Infinity V500, Evita V600 and Evita V800 do not send this code.

Code	Code page	Alarm description	Alarm phrase	New with version	Used by device	Remarks
E8H	1	Tidal volume > high limit	TIDAL VOL HI	01.00	Perseus A500, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Savina 300, Zeus IE	
E9H	1	O2 alarm disabled	O2 ALARM OFF	01.00	Savina 300	
EBH	1	CO2 measurement values can temporarily not be measured with accuracy complying to ISO standard.	CO2 INACC.	01.07	Perseus A500, Atlan, Evita V600, Evita V800, Babylog VN600, Babylog VN800, Primus Family, Zeus IE (added with SW 2.00 for Perseus A500, added with SW 2.00 for Evita V800, Evita V600, Babylog VN800 and Babylog VN600, added with SW 2.00 for Zeus IE, added with SW 4.52 for the Primus Family)	
ECH	1	Gas supply failure	GAS FAILURE	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800	
EDH	1	Check N2O supply	N2O SUPPLY ?	01.00	Perseus A500, Atlan, Primus Family, Zeus IE	
EEH	1	Two agents detected	2nd AGENT	01.00	Perseus A500, Atlan, Primus Family, Zeus IE	

Code	Code page	Alarm description	Alarm phrase	New with version	Used by device	Remarks
EFH	1	Power fail	POWER FAIL	01.00	Perseus A500, Atlan, Primus Family, Zeus IE, Carina, Fabius Family	
F0H	1	Skin temperature difference high	DTEMP HIGH	01.08	Babyleo TN500, Isolette 8000 plus, Babyroo TN300	
F1H	1	Check setting of Pmax	P MAX?	01.01	Primus Family, Carina	
F2H	1	Internal system fault	SYSTEM FAULT	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Carina, Babyleo TN500, Isolette 8000 plus, Babyroo TN300	
F3H	1	No fresh gas	NO FRESHGAS	01.00	Perseus A500, Atlan, Apollo, Zeus IE, Fabius Family except Fabius plus, Fabius Plus XL	Edition 18: From the Primus Family only the Apollo sends this code.
F4H	1	O2 safety flow activated	SAFETY O2 ON	01.00	Perseus A500, Atlan, Primus Family, Zeus IE	
F5H	1	Inspiratory MAC > high limit	INSP MAC HI	01.00	Perseus A500, Atlan, Zeus IE	
F6H	1	Skin temperature difference low	DTEMP LOW	01.08	Babyleo TN500, Isolette 8000 plus	
F7H	1	Inspiratory CO2 alarm disabled	FICO2 '@ OFF	01.01	Primus Family	

Code	Code page	Alarm description	Alarm phrase	New with version	Used by device	Remarks
F8H	1	Continuous high pressure	CONT PRES	01.00	Perseus A500, Atlan, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300, Carina, Fabius Family	
FAH	1	Pressure threshold low	THRESHOLD LO	01.05	Fabius Family	
FBH	1	Check APL valve	APL VALVE ?	01.05	Fabius Family	
FCH	1	CO2 not calibrated	CO2 NOT CAL	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Oxylog 3000 plus, Oxylog VE300	
FDH	1	Battery inoperable	BATTERY ERR	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Zeus IE, Oxylog 3000 plus, Oxylog VE300, Fabius Family (added with SW 2.00 for Perseus A500)	
FEH	1	Check cooling	COOLING ?	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300	

Code	Code page	Alarm description	Alarm phrase	New with version	Used by device	Remarks
FFH	1	Disconnection ventilator	DISCONNECT	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800	
01H	2	NIBP did not measure	NBP NO RSLT	01.01	Zeus IE	
03H	2	NIBP max. inflation press. exceeded	OVERPRESS	01.01	Zeus IE	
07H	2	NIBP measurement max time exceeded	NBP TIMEOUT	01.01	Zeus IE	
0BH	2	NIBP tube blocked	NBP TUBE ?	01.01	Zeus IE	
0CH	2	Humidifier inoperable	HUM INOP	01.08	Babyleo TN500, Isolette 8000 plus	
0DH	2	Humidity low	HUMIDITY LOW	01.08	Babyleo TN500, Isolette 8000 plus	
10H	2	N2O supply pressure > high limit	N2O PRESS HI	01.00	Perseus A500, Atlan	
11H	2	Volatile agent dosage error	VA ERR	01.01	Zeus IE	
12H	2	Volatile agent dosage error and mixer error	VA+MIX ERR	01.01	Zeus IE	
14H	2	CPAP pressure not established due to leak. CPAP setting has been disabled.	CPAP DEACT.	01.03	Perseus A500 (added with SW 1.11 for Perseus A500)	
15H	2	Central skin temperature < low limit	CTEMP LOW	01.08	Babyleo TN500, Isolette 8000 plus, Babyroo TN300	
16H	2	Central skin temperature > high limit	CTEMP HIGH	01.08	Babyleo TN500, Isolette 8000 plus, Babyroo TN300	
17H	2	Central skin temperature failure	CTEMP ERR	01.08	Babyleo TN500, Isolette 8000 plus, Babyroo TN300	
18H	2	Peripheral skin temperature < low limit	PTEMP LOW	01.08	Babyleo TN500, Isolette 8000 plus, Babyroo TN300	
19H	2	Peripheral skin temperature > high limit	PTEMP HIGH	01.08	Babyleo TN500, Isolette 8000 plus, Babyroo TN300	
1AH	2	Peripheral skin temperature failure	PTEMP ERR	01.08	Babyleo TN500, Isolette 8000 plus, Babyroo TN300	
1BH	2	Ambient temperature measurement failure	AMB TEMP ERR	01.09	Babyleo TN500, Babyroo TN300	
31H	2	O2 cylinder pressure low	O2CYLIND.LOW	01.01	Primus Family	

Code	Code page	Alarm description	Alarm phrase	New with version	Used by device	Remarks
32H	2	O2 cylinder empty or cylinder pressure low	O2 CYL. ?	01.00	Perseus A500, Atlan, Primus Family, Zeus IE	
33H	2	O2 cylinder not connected	O2CYL.DISCON	01.01	not used	Edition 18: The Primus Family does not send this code.
34H	2	N2O cylinder empty or cylinder pressure low	N2O CYL.?	01.00	Perseus A500, Atlan, Primus Family, Zeus IE	
35H	2	N2O delivery failure	NO N2O	01.00	Perseus A500, Atlan, Apollo, Zeus IE	Edition 6: Primus and Primus IE do not send this code.
36H	2	O2 delivery failure	NO OXYGEN	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Apollo, Zeus IE, Carina, Babyleo TN500	Edition 6: Primus and Primus IE do not send this code.
37H	2	AIR delivery failure	NO AIR	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Apollo, Zeus IE	Edition 6: Primus and Primus IE do not send this code.
38H	2	Set fresh-gas flow not attained	FG LIMITED	01.00	Perseus A500, Atlan, Primus Family	
39H	2	Internal/external switch valve inoperable	FG EXTERN?	01.00	Perseus A500, Atlan, Primus Family	

Code	Code page	Alarm description	Alarm phrase	New with version	Used by device	Remarks
3AH	2	Insp. N2O > high limit	% N2O HIGH	01.00	Perseus A500, Atlan, Primus Family, Zeus IE	
3BH	2	Ambient pressure measurement disturbed or inoperable	AMB PRESS ?	01.00	Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800 (added with SW 2.10 for Atlan)	
3CH	2	SDC certificate cannot be used	SDC CERT ERR	01.15	Evita V600, Evita V800, Babylog VN600, Babylog VN800 (added with SW 3.00 for Evita V800, Evita V600, Babylog VN800 and Babylog VN600)	
40H	2	Agent dosage inoperable	VOLAT ERR	01.01	Zeus IE	
41H	2	Agent reservoir low	VOLAT SUPPLY	01.00	Perseus A500, Zeus IE	
42H	2	CO2 sample gas line disconnected	CO2-LINE ?	01.00	Perseus A500, Atlan, Zeus IE	
43H	2	Check setting flow	CHECK FLOW	01.01	Oxylog 3000 plus	
44H	2	Check setting FiO2	CHECK FIO2	01.01	Oxylog 3000 plus	
45H	2	High-flow - switch to fresh gas control mode!	FG-FLOW HIGH	01.01	Zeus IE	
4FH	2	Network communication lost	NETWORK ERR	01.15	Evita V600, Evita V800, Babylog VN600, Babylog VN800 (added with SW 3.00 for Evita V800, Evita V600, Babylog VN800 and Babylog VN600)	
50H	2	Oxygen module inoperable	O2 MOD INOP	01.09	Babyleo TN500	
59H	2	No internal battery detected	INT.BATTERY?	01.01	Oxylog 3000 plus	

Code	Code page	Alarm description	Alarm phrase	New with version	Used by device	Remarks
5AH	2	Power supply by battery	BATTERY ON	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800	
5BH	2	No nebuliser	NO NEBULIS.	01.00	Savina 300	
5CH	2	Battery less than 2 min left	BATT. < 2MIN	01.00	Evita V300, Evita Infinity V500, Babylog VN500, Oxylog 3000 plus, Oxylog VE300, Carina (removed with SW 2.40 for Babylog VN500, removed with SW 2.40 for Evita Infinity V500 and Evita V300)	Edition 7: beginning with SW 2.40 Evita Infinity V500, Evita V300 and Babylog VN500 no longer use this code
5DH	2	Battery less than 5 min left	BATT. < 5MIN	01.04	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800 (added with SW 2.40 for Babylog VN500, added with SW 2.40 for Evita Infinity V500 and Evita V300)	
66H	2	Apnea alarm disabled	APN ALRM OFF	01.00	Savina 300, Zeus IE, Fabius Family	
67H	2	Minute volume low alarm disabled	MV LOW OFF	01.00	Savina 300, Primus Family	
68H	2	Tidal volume high alarm disabled	VT HIGH OFF	01.00	Savina 300	

Code	Code page	Alarm description	Alarm phrase	New with version	Used by device	Remarks
69H	2	Distribution of alarms via SDC not possible	DAS ERR	01.15	Evita V600, Evita V800, Babylog VN600, Babylog VN800 (added with SW 3.00 for Evita V800, Evita V600, Babylog VN800 and Babylog VN600)	
6AH	2	Tube obstructed	TUBE OBSTRUC	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Oxylog VE300	
6CH	2	Check microfilter	MICROFILTER?	01.00	Savina 300, Carina	
6DH	2	Power supply by external DC	EXTERN DC ON	01.00	Savina 300	
6EH	2	Ambient pressure high	AMB PRESS HI	01.00	Savina 300	
6FH	2	Ambient pressure low	AMB PRESS LO	01.00	Savina 300, Oxylog VE300	
74H	2	Right atrial pressure < low limit	RA LO	01.01	Zeus IE	
75H	2	Right atrial pressure > high limit	RA HI	01.01	Zeus IE	
78H	2	Problems with right atrial pressure measurement	RA ERR	01.01	Zeus IE	
7CH	2	Mattress temperature high	MAT TEMP HI	01.09	Babyleo TN500, Babyroo TN300	
7DH	2	Mattress temperature sensor inoperable	MAT TEMP ERR	01.13	Babyroo TN300	
7FH	2	Check patient	CHK PATIENT	01.09	Babyleo TN500, Babyroo TN300	
80H	2	Radiant warmer faulty	RAD WRM INOP	01.13	Babyroo TN300	
81H	2	Mattress temperature low	MAT TEMP LO	01.13	Babyroo TN300	
8DH	2	Breathing system disconnected	VENT ASSEMBL	01.00	Perseus A500, Atlan, Primus Family	
90H	2	Check neonatal flow sensor	NEO FLOW?	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800	

Code	Code page	Alarm description	Alarm phrase	New with version	Used by device	Remarks
91H	2	Loss of data	LOSS OF DATA	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Primus Family, Oxylog 3000 plus, Oxylog VE300, Babyleo TN500, Isolette 8000 plus, Babyroo TN300	
93H	2	Apnea ventilation	APNEA VENT	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Oxylog 3000 plus, Oxylog VE300, Carina, Fabius Family (added with SW 2.00 for Oxylog VE300)	
94H	2	Check ventilator	CHECK VENT	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300	
95H	2	Ventilator standby	VENT STANDBY	01.00	Savina 300	
96H	2	Problems with PEEP control	PEEP ERR	01.00	Savina 300	
98H	2	Nebulizer is activated	NEBULIZER ON	01.00	Savina 300	

Code	Code page	Alarm description	Alarm phrase	New with version	Used by device	Remarks
99H	2	Inspiration hold aborted	INSPHOLD END	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300	
9CH	2	Leakage detected	LEAKAGE	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Oxylog 3000 plus, Oxylog VE300, Carina	
9DH	2	Backup ventilation	BACKUP VENT	01.00	Perseus A500, Atlan	
9EH	2	Expiration hold aborted	EXPHOLD END	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Savina 300	
9FH	2	Nebulization terminated	NEBULIZ. OFF	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800	
A0H	2	Ventilator not in locked Position	VENT.UNLOCKD	01.01	Primus Family	

Code	Code page	Alarm description	Alarm phrase	New with version	Used by device	Remarks
A1H	2	Problems with power supply	PWR SPLY ERR	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Primus Family, Oxylog 3000 plus, Oxylog VE300, Carina	
A2H	2	Expiratory halothan > high limit	EXP.HAL HIGH	01.01	not used	Edition 18: The Primus Family does not send this code.
A3H	2	Expiratory enflurane > high limit	EXP.ENF HIGH	01.01	not used	Edition 18: The Primus Family does not send this code.
A4H	2	Expiratory isoflurane > high limit	EXP.ISO HIGH	01.01	not used	Edition 18: The Primus Family does not send this code.
A5H	2	Expiratory desflurane > high limit	EXP.DES HIGH	01.01	not used	Edition 18: The Primus Family does not send this code.
A6H	2	Expiratory sevoflurane > high limit	EXP.SEV HIGH	01.01	not used	Edition 18: The Primus Family does not send this code.
A7H	2	Set tidal volume not attained	TIDAL VOL.?	01.00	Perseus A500, Atlan, Primus Family	
A8H	2	Inspiratory flow measurement inoperable	INSP VOL ERR	01.00	Perseus A500, Atlan, Primus Family, Oxylog 3000 plus	
A9H	2	Setting could not be performed	SET.CANCELED	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Primus Family, Oxylog 3000 plus, Oxylog VE300, Carina	
ABH	2	Confirm settings	CONFIRM SET.	01.01	Oxylog 3000 plus	

Code	Code page	Alarm description	Alarm phrase	New with version	Used by device	Remarks
ACH	2	Fresh-gas flow too high	FG TOO HIGH	01.00	Perseus A500, Atlan, Apollo	Edition 18: From the Primus Family only the Apollo sends this code.
ADH	2	Fresh-gas flow active	FG ACTIVE	01.00	Perseus A500, Atlan, Apollo	Edition 18: From the Primus Family only the Apollo sends this code.
AEH	2	Air trapping	AIR TRAPPING	01.01	Zeus IE	
AFH	2	Oxygen cylinder open	O2 CYL OPEN	01.00	Apollo	Edition 6: Perseus A500, Primus and Primus IE do not send this code.
B0H	2	N2O cylinder open	N2O CYL OPEN	01.00	Apollo	Edition 6: Perseus A500, Primus and Primus IE do not send this code.
B1H	2	Air cylinder open	AIR CYL OPEN	01.00	Apollo	Edition 6: Perseus A500, Primus and Primus IE do not send this code.
B2H	2	N2O cylinder sensor not connected	N2OCYL.SENS?	01.00	Perseus A500, Atlan, Primus Family	
B3H	2	Air cylinder sensor not connected	AIRCYL.SENS?	01.00	Perseus A500, Atlan, Primus Family	
B4H	2	O2 cylinder sensor not connected	O2 CYL.SENS?	01.00	Perseus A500, Atlan, Primus Family	
B5H	2	Air cylinder empty or cylinder pressure low	AIR CYL.?	01.00	Perseus A500, Atlan, Primus Family, Zeus IE	
B6H	2	Inspiratory resistance >> resistance of test	IN.RES. HIGH	01.01	Zeus IE	
B7H	2	Expiratory resistance << resistance of test	EXP.RES.HIGH	01.01	Zeus IE	

Code	Code page	Alarm description	Alarm phrase	New with version	Used by device	Remarks
B8H	2	A setting, alarm limit or ventilation mode was changed, but the change was not confirmed. If necessary, the user can adjust and confirm the setting again with the rotary knob.	NO CONFIRM.	01.05	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Babyleo TN500, Babyroo TN300 (added with SW 2.42 for Evita Infinity V500, Evita V300 and Babylog VN500, added with SW 4.20 for Savina 300)	
B9H	2	Automatic uptake control not operable	PBAG INOP	01.01	Zeus IE	
BAH	2	Delivered volume greater set tidal volume due to minimum required PIP	PMIN REACHED	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Zeus IE (added with SW 2.00 for Atlan, removed with SW 3.00 for Evita V800 and Evita V600)	
BBH	2	Insp. N2O < low limit	% N2O LOW	01.01	Zeus IE	
BCH	2	Detected agent differs from selected type	WRONG AGENT	01.01	Zeus IE	
BDH	2	Insp. vol. agent > high limit	INSP AGT.HI.	01.01	Zeus IE	
BEH	2	DIVA module -I- or -II- is warming up	DIVA WARMUP	01.01	Zeus IE	
BFH	2	Relearning ECG pattern	ECG RELEARN	01.01	Zeus IE	
C0H	2	Unable to relearn ECG pattern	ECG RLRN.ERR	01.01	Zeus IE	
C1H	2	ECG inoperable	ECG INOP	01.01	Zeus IE	
C2H	2	RL-electrode failure	ECG ELEC. RL	01.01	Zeus IE	
C3H	2	V-plus electrode failure	ECG ELEC. V+	01.01	Zeus IE	
C4H	2	Cardiac output blood temperature > high limit	BLOOD TEMP"#	01.01	Zeus IE	
C5H	2	Cardiac output blood temperature > low limit	BLOOD TEMP\$&	01.01	Zeus IE	
C6H	2	All alarms disabled	ALL ALRM SUS	01.01	Zeus IE	

Code	Code page	Alarm description	Alarm phrase	New with version	Used by device	Remarks
C7H	2	Air fresh-gas flow measurement inoperable	FG AIR SENS?	01.01	Perseus A500, Atlan, Apollo (added with SW 1.11 for Perseus A500)	Edition 7: Primus and Primus IE do not send this code.
C8H	2	O2 fresh-gas flow measurement inoperable	FG O2 SENS?	01.01	Perseus A500, Atlan, Apollo (added with SW 1.11 for Perseus A500)	Edition 7: Primus and Primus IE do not send this code.
C9H	2	N2O fresh-gas flow measurement inoperable	FG N2O SENS?	01.01	Perseus A500, Atlan, Apollo (added with SW 1.11 for Perseus A500)	Edition 7: Primus and Primus IE do not send this code.
CAH	2	No air supply	NO AIR	01.00	Babylog VN500, Babylog VN600, Babylog VN800	
CBH	2	No N2O supply (deprecated, will not be used by MEDIBUS.X devices)	NO N2O	01.00	not used	
CCH	2	PEEP < low limit	PEEP LOW	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800	
CDH	2	Set inspiratory pressure not attained	PINSP ?	01.00	Perseus A500, Atlan, Oxylog VE300, Fabius Family	
CEH	2	Hose kinked	HOSE KINKED	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800	
CFH	2	Set tidal volume too high for connected hose	VT HIGH HOSE	01.01	Oxylog 3000 plus	
D0H	2	Pressure relief valve opened	PRESS RELIEF	01.00	Perseus A500, Atlan, Primus Family	
D1H	2	Hose system defect	HOSE ERROR	01.00	Evita V300, Evita Infinity V500, Babylog VN500	

Code	Code page	Alarm description	Alarm phrase	New with version	Used by device	Remarks
D2H	2	Check MAC value	MAC LOW?	01.00	Perseus A500, Atlan, Primus Family	
D3H	2	Water trap expired?	WATERTR.OLD?	01.00	Perseus A500, Atlan, Primus Infinity Empowered	Edition 6: Primus and Apollo do not send this code.
D4H	2	CO2 absorber depleted?	ABSORB. OLD?	01.00	Perseus A500, Atlan, Primus Infinity Empowered	Edition 6: Primus and Apollo do not send this code.
D5H	2	CO2 absorber present?	ABS.PRESENT?	01.00	Perseus A500, Atlan, Primus Infinity Empowered	Edition 6: Primus and Apollo do not send this code.
D6H	2	Hoses interchanged?	HOSES MIXED?	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Babylog VN500, Primus Infinity Empowered	Edition 6: Primus and Apollo do not send this code.
D7H	2	Hoses incompatible?	WRONG HOSES?	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Babylog VN500, Primus Infinity Empowered, Oxylog 3000 plus, Oxylog VE300	Edition 6: Primus and Apollo do not send this code.
D8H	2	Accessory ID detection functions inoperable	ID-FUNC-INOP	01.00	Evita V300, Evita Infinity V500, Babylog VN500, Primus Infinity Empowered	Edition 19: Perseus A500 and Atlan do not use this code. Edition 6: Primus and Apollo do not send this code.
D9H	2	Breathing circuit expired?	HOSE OLD?	01.00	Perseus A500, Atlan, Primus Infinity Empowered	Edition 6: Primus and Apollo do not send this code.
DAH	2	Set expiration time not attainable	EXP TIME ERR	01.00	Evita V300, Evita Infinity V500, Babylog VN500, Oxylog 3000 plus	
DBH	2	NMT TOF > high limit	NMT TOF.HIGH	01.01	Zeus IE	

Code	Code page	Alarm description	Alarm phrase	New with version	Used by device	Remarks
DCH	2	NMT TOF < low limit	NMT TOF.LOW	01.01	Zeus IE	
DDH	2	NMT initializing/NMT single out of range	NMT INOP	01.01	Zeus IE	
DEH	2	NMT check sensors and/or electrodes	NMT SENS?	01.01	Zeus IE	
DFH	2	NMT temperature out of range	NMT TEMP?	01.01	Zeus IE	
E0H	2	Bispectral index high	BIS.HIGH	01.01	Zeus IE	
E1H	2	Bispectral index low	BIS.LOW	01.01	Zeus IE	
E2H	2	Bispectral calculation not possible	BIS INOP	01.01	Zeus IE	
E3H	2	Check BIS sensor	BIS SENS?	01.01	Zeus IE	
E4H	2	BISX out of range	BIS?	01.01	Zeus IE	
E5H	2	SmartCare: patient session aborted	SC ABORTED	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Zeus IE (added with SW 1.10 for Zeus IE)	
E6H	2	SmartCare: inoperable / patient session terminated	SC INOP	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800	
E7H	2	SmartCare: Central Hypoventilation	CENTRAL HYPO	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800	
E8H	2	SmartCare: Persistent Tachypnea	PERS TACHYP	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800	
E9H	2	SmartCare: Unexplained Hyperventilation	UNEXPL HYPER	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800	
EAH	2	SmartCare: Reduce PEEP if possible	PEEP REDUCIB	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800	
EBH	2	SmartCare: SBT successful	CONS SEPARAT	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800	

Code	Code page	Alarm description	Alarm phrase	New with version	Used by device	Remarks
ECH	2	Ejector inoperable	EJECTOR INOP	01.00	Babylog VN500, Babylog VN600, Babylog VN800	
EDH	2	Constant CO2 offset	CONSTANT CO2	01.10	Oxylog VE300	
EFH	2	Autoflow VT too low	AF VT LOW	01.01	Oxylog 3000 plus	
F1H	2	Pilot lines switched	PILOTLINE SW	01.01	Oxylog 3000 plus, Oxylog VE300	
F2H	2	Breathing hoses missing	HOSE MISSING	01.01	Primus Infinity Empowered	Edition 6: Primus and Apollo do not send this code.
F3H	2	SmartCare: Reduce FiO2 if possible	FIO2 REDUCIB	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800	
F4H	2	Mean Airway Pressure < low limit	MAP LOW	01.00	Babylog VN500, Babylog VN600, Babylog VN800	
F7H	2	Airway pressure high (HFO)	PAW HIGH HF	01.00	Babylog VN500, Babylog VN600, Babylog VN800	
F8H	2	Tidal volume < low Limit	TIDAL VOL LO	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300	
F9H	2	Plow > high limit	PLOW HIGH	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800	
FAH	2	Plow < low limit	PLOW LOW	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800	

Code	Code page	Alarm description	Alarm phrase	New with version	Used by device	Remarks
FBH	2	Delivered volume limited by VTmax	VT LIMITED	01.00	Babylog VN500, Babylog VN600, Babylog VN800 (removed with SW 3.00 for Babylog VN800 and Babylog VN600)	
FCH	2	Air entrainment	AIR ENTRAIN	01.00	Perseus A500, Atlan	
FDH	2	Ventilation pause exceeded	VENT PAUSE?	01.00	Perseus A500, Atlan	
FEH	2	Power supply by internal battery	INT BATT ON	01.00	Savina 300, Oxylog 3000 plus, Carina	
FFH	2	In O2 therapy the set flow cannot be reached due to high resistance.	FLOW LOW	01.08	Evita V600, Evita V800, Babylog VN600, Babylog VN800, Savina 300 (added with SW 2.00 for Evita V800, Evita V600, Babylog VN800 and Babylog VN600, added with SW 5.00 for Savina 300)	

Text messages

Code	Description	Text message (Note: '_' represents the space character 20H)	New with version	Used by device	Remarks
01H	Mode VC-CMV (volume-controlled, continuous mandatory ventilation)	_Mode VC-CMV	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300, Fabius Family	
02H	Mode VC-AC (volume-controlled, assist controlled ventilation)	_Mode VC-AC	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Savina 300, Oxylog 3000 plus, Oxylog VE300, Carina	
06H	Mode VC-SIMV (volume-controlled, synchronized intermittent mandatory ventilation)	_Mode VC-SIMV	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300, Carina, Fabius Family	

Code	Description	Text message (Note: ' ' represents the space character 20H)	New with version	Used by device	Remarks
0AH	Mode SPN-CPAP (spontaneous ventilation, continuous positive airway pressure)	_Mode SPN-CPAP	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300, Carina, Fabius Family	
0CH	Mode VC-MMV (volume-controlled, mandatory minute ventilation)	_Mode VC-MMV	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Savina 300	
0DH	Mode PC-MMV (pressure-controlled, mandatory minute ventilation)	_Mode PC-MMV	01.00	Babylog VN500, Babylog VN600, Babylog VN800	
0EH	Mode PC-BIPAP (pressure-controlled, bi-phasic positive airway pressure)	_Mode PC-BIPAP	01.00	Perseus A500, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Savina 300, Zeus IE, Oxylog 3000 plus, Oxylog VE300, Carina (added with SW 1.10 for Oxylog VE300)	

Code	Description	Text message (Note: ' ' represents the space character 20H)	New with version	Used by device	Remarks
11H	Automatic apnea ventilation active	APNEA VENTILATION	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Oxylog 3000 plus, Oxylog VE300, Carina (added with SW 2.00 for Oxylog VE300)	
12H	Device in PRODUCT TEST mode	Mode DS	01.00	Perseus A500, Atlan, Babyleo TN500, Isolette 8000 plus, Babyroo TN300	Edition 6: Evita Infinity V500 and Babylog VN500 do not send this code.
13H	Air mode is active	Air Mode	01.08	Babyleo TN500, Isolette 8000 plus	
14H	Skin mode is active	Skin Mode	01.08	Babyleo TN500, Isolette 8000 plus, Babyroo TN300	
15H	Oxygen control is active	Oxygen Control active	01.08	Babyleo TN500, Isolette 8000 plus	
16H	Humidity control is active	Humidity Control active	01.08	Babyleo TN500, Isolette 8000 plus	
18H	Mode PC-SIMV (pressure-controlled, synchronized intermittent mandatory ventilation)	_Mode PC-SIMV	01.00	Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Primus Family	
1AH	Mode PC-APRV (pressure-controlled, airway pressure release ventilation)	_Mode PC-APRV	01.00	Perseus A500, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300	

Code	Description	Text message (Note: ' ' represents the space character 20H)	New with version	Used by device	Remarks
1BH	Mode PC-HFO (pressure-controlled, high frequency oscillation)	_Mode PC-HFO	01.00	Babylog VN500, Babylog VN600, Babylog VN800	
1CH	Suction maneuver active	Suction maneuver active	01.01	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300	
1DH	Low-flow maneuver active	Low-flow maneuver active	01.01	Evita Infinity V500, Evita V600, Evita V800	
1EH	Ventilator is in standby mode	Ventilator STANDBY	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog VE300, Carina	
20H	Device is in adult mode	Mode Adults	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Zeus IE, Oxylog VE300 (added with SW 2.02 for Zeus IE)	

Code	Description	Text message (Note: ' ' represents the space character 20H)	New with version	Used by device	Remarks
21H	Device is in neonatal mode	Mode Neonates	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Zeus IE (added with SW 2.02 for Zeus IE)	
22H	Selected CO2 unit is mmHg	mmHg	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300 (added with SW 4.20 for Savina 300)	
23H	Selected CO2 unit is kPa	kPa	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300 (added with SW 4.20 for Savina 300)	

Code	Description	Text message (Note: ' ' represents the space character 20H)	New with version	Used by device	Remarks
24H	Selected CO2 unit is %	%	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300 (added with SW 4.20 for Savina 300)	
25H	Halothane selected or detected	Anesthesia gas HALOTHANE	01.00	Perseus A500, Atlan, Primus Family, Zeus IE, Fabius Family	
26H	Enflurane selected or detected	Anesthesia gas ENFLURANE	01.00	Perseus A500, Atlan, Primus Family, Zeus IE, Fabius Family	
27H	Isoflurane selected or detected	Anesthesia gas ISOFLURANE	01.00	Perseus A500, Atlan, Primus Family, Zeus IE, Fabius Family	
28H	Desflurane selected or detected	Anesthesia gas DESFLURANE	01.00	Perseus A500, Atlan, Primus Family, Zeus IE, Fabius Family	
29H	Sevoflurane selected or detected	Anesthesia gas SEVOFLURANE	01.00	Perseus A500, Atlan, Primus Family, Zeus IE, Fabius Family	
2AH	No anesthesia agent selected nor detected	No anesthesia gas	01.00	Perseus A500, Atlan, Primus Family, Zeus IE	

Code	Description	Text message (Note: ' ' represents the space character 20H)	New with version	Used by device	Remarks
2BH	Ventilation mode man./spont.	_Mode MAN/SPONT	01.00	Perseus A500, Atlan, Primus Family, Zeus IE, Fabius Family	
2CH	Selected language (variable text string transmitted with this code, see "MEDIBUS.X, Rules and Standards for Implementation" for further information)		01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Fabius Family	
35H	Mode SPN-PPS (spontaneous, proportional pressure support)	_Mode SPN-PPS	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800	
36H	Ventilation mode FRESH GAS EXTERNAL	_Mode FRESH GAS EXT.	01.00	Perseus A500, Atlan, Primus Family, Zeus IE	
37H	Selected carrier gas is air	Carrier Gas AIR	01.00	Perseus A500, Atlan, Primus, Primus Infinity Empowered, Zeus IE	Edition 6: Apollo does not send this code.
38H	Selected carrier gas is N2O	Carrier Gas N2O	01.00	Perseus A500, Atlan, Primus, Primus Infinity Empowered, Zeus IE	Edition 6: Apollo does not send this code.

Code	Description	Text message (Note: ' ' represents the space character 20H)	New with version	Used by device	Remarks
3AH	Device is in pediatric mode	Mode Pediatrics	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Zeus IE, Oxylog VE300 (added with SW 2.02 for Zeus IE)	
3BH	Manual mode is active	Manual Mode	01.09	Babyleo TN500, Babyroo TN300	
3CH	Open care therapy is active	Open Care Therapy	01.09	Babyleo TN500	
3DH	Photo therapy is active	Photo Therapy on	01.10	BiliLux	
3EH	Mode PC-PSV (pressure controlled, pressure support ventilation)	_Mode PC-PSV	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800	
42H	Patient is a female	female	01.00	Perseus A500, Zeus IE (removed with SW 1.11 for Perseus A500)	Edition 7: beginning with SW 1.11 Perseus A500 is no longer using this code
43H	Patient is a male	male	01.00	Perseus A500, Zeus IE (removed with SW 1.11 for Perseus A500)	Edition 7: beginning with SW 1.11 Perseus A500 is no longer using this code
45H	End-expiratory control of anesthesia gas active	End-exp. controlled Agent	01.01	Zeus IE	
47H	Mode PC-AC (pressure controlled, assisted-controlled)	_Mode PC-AC	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Carina	

Code	Description	Text message (Note: ' ' represents the space character 20H)	New with version	Used by device	Remarks
48H	Device configured for intubated patient ventilation	IV - Invasive Ventilation	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Oxylog 3000 plus, Oxylog VE300, Carina	
49H	Device configured for mask ventilation	NIV - Non-Invasive Ventilation	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Oxylog 3000 plus, Oxylog VE300, Carina	
4AH	2nd agent halothane selected or detected	2nd Anesthesia gas HALOTHANE	01.00	Perseus A500, Atlan, Primus Family, Zeus IE, Fabius Family	
4BH	2nd agent enflurane selected or detected	2nd Anesthesia gas ENFLURANE	01.00	Perseus A500, Atlan, Primus Family, Zeus IE, Fabius Family	
4CH	2nd agent isoflurane selected or detected	2nd Anesthesia gas ISOFLURANE	01.00	Perseus A500, Atlan, Primus Family, Zeus IE, Fabius Family	
4DH	2nd agent desflurane selected or detected	2nd Anesthesia gas DESFLURANE	01.00	Perseus A500, Atlan, Primus Family, Zeus IE, Fabius Family	

Code	Description	Text message (Note: ' ' represents the space character 20H)	New with version	Used by device	Remarks
4EH	2nd agent sevoflurane selected or detected	2nd Anesthesia gas SEVOFLURANE	01.00	Perseus A500, Atlan, Primus Family, Zeus IE, Fabius Family	
4FH	No 2nd anesthesia agent selected nor detected	No 2nd anesthesia gas	01.00	Perseus A500, Atlan, Primus Family, Zeus IE	
50H	Mattress temperature control is active	Mattress Temp Control active	01.09	Babyleo TN500, Babyroo TN300	
51H	Kangaroo mode is active	Kangaroo Mode	01.08	Babyleo TN500, Isolette 8000 plus, Babyroo TN300	
53H	Device is performing leakage test	performing LEAKAGE TEST	01.01	Primus Family, Zeus IE	
54H	Device is in standby mode	Mode STANDBY	01.00	Perseus A500, Atlan, Primus Family, Zeus IE, Fabius Family	
55H	Low pressure oxygenation active	LPO - Low Pressure Oxygenation	01.00	Savina 300	
56H	Selected agent unit is kPa	kPa	01.00	Perseus A500, Atlan, Primus Family	
57H	Selected agent unit is %	%	01.00	Perseus A500, Atlan, Primus Family	
58H	HLM mode active	HLM Mode active	01.00	Perseus A500, Atlan, Primus Family, Zeus IE	

Code	Description	Text message (Note: ' ' represents the space character 20H)	New with version	Used by device	Remarks
5CH	Pressure support added to current ventilation mode	/PS	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Oxylog VE300, Carina, Fabius Family	
5DH	Audio pause active	Audio pause active	01.01	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300 (added with SW 2.00 for Atlan, added with SW 2.03 for Perseus A500)	Edition 20: Former edition states that this code is used by SW 1.02 for Atlan. This code is not be used by SW 1.02 for Atlan. It has been added for Atlan with SW 2.00
5EH	AutoFlow added to current ventilation mode	/AF	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Savina 300, Primus Family, Zeus IE, Oxylog 3000 plus, Carina	
5FH	Fresh-gas control of anaesthesia gas active	Fresh Gas Control	01.01	Zeus IE	

Code	Description	Text message (Note: ' ' represents the space character 20H)	New with version	Used by device	Remarks
60H	First name of patient (variable text string transmitted with this code, see "MEDIBUS.X, Rules and Standards for Implementation" for further information)		01.01	Zeus IE	
61H	Last name of patient (variable text string transmitted with this code, see "MEDIBUS.X, Rules and Standards for Implementation" for further information)		01.01	Zeus IE	
62H	Patient monitoring is in simulation mode	SIMULATION ON	01.01	Zeus IE	
63H	Ventilator was switched off by emergence switch "Anesthesia ventilator OFF"	Therapy off	01.01	Zeus IE	
64H	Trigger level normal	Trigger level normal	01.01	Carina	
65H	Trigger level sensitive	Trigger level sensitive	01.01	Carina	
67H	Apnea detection disabled	Apnea detection disabled	01.01	Carina	
68H	Tube type endotracheal	Endotracheal	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300 (added with SW 5.00 for Savina 300)	
69H	Volume guarantee added to current ventilation mode	/VG	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Carina	
6AH	Mode PC-CMV (pressure controlled, continuous mandatory ventilation)	_Mode PC-CMV	01.00	Perseus A500, Atlan, Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Primus Family, Zeus IE, Fabius Family	

Code	Description	Text message (Note: ' ' represents the space character 20H)	New with version	Used by device	Remarks
6BH	Mode O2 therapy	_Mode O2 Therapy	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300 (added with SW 5.00 for Savina 300)	
6CH	Abdominal external trigger	Abdominal External Trigger	01.00		
6DH	Automatic tube compensation (ATC) added to current ventilation mode	/ATC	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300 (added with SW 5.00 for Savina 300)	
6EH	Expiratory automatic tube compensation (ATC) enabled	Expiratory ATC	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800	
6FH	Inspiratory automatic tube compensation (ATC) enabled	Inspiratory ATC	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800	
70H	Tube type tracheostoma	Tracheostoma	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Savina 300 (added with SW 5.00 for Savina 300)	

Code	Description	Text message (Note: ' ' represents the space character 20H)	New with version	Used by device	Remarks
71H	Apnea ventilation enabled	Apnea Ventilation enabled	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Oxylog 3000 plus, Oxylog VE300 (added with SW 2.00 for Oxylog VE300)	
72H	Apnea ventilation auto-return enabled	Apnea Ventilation Auto Return	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800	
73H	Screen is locked	Screen locked	01.01	Carina	
74H	Volume support added to current ventilation mode	/VS	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800	
75H	AutoRelease added to current ventilation mode	/AutoRelease	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800	
76H	Active humid unheated	Active humid unheated	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800 (removed with SW 3.00 for Evita V800, Evita V600, Babylog VN800 and Babylog VN600)	Edition 22: Starting with SW 3.00 Evita V600, Evita V800, Babylog VN600 and Babylog VN800 use code AAH instead.

Code	Description	Text message (Note: ' ' represents the space character 20H)	New with version	Used by device	Remarks
77H	Active humid heated	Active humid heated	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800 (removed with SW 3.00 for Evita V800, Evita V600, Babylog VN800 and Babylog VN600)	Edition 22: Starting with SW 3.00 Evita V600, Evita V800, Babylog VN600 and Babylog VN800 use code ABH instead.
78H	HME / Filter	HME / Filter	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800, Oxylog 3000 plus, Oxylog VE300	
79H	SmartCare therapy phase: adapting (Note: only one SmartCare therapy phase can be active at a time)	SC Therapy phase: Adapting	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800	
7AH	SmartCare therapy phase: observing	SC Therapy phase: Observing	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800	
7BH	SmartCare therapy phase: maintaining	SC Therapy phase: Maintaining	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800	
7CH	SmartCare diagnosis: normal ventilation (Note: only one SmartCare diagnosis can be active at a time)	Diagnosis: Normal Ventilation	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800	
7DH	SmartCare diagnosis: hyperventilation	Diagnosis: Hyperventilation	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800	
7EH	SmartCare diagnosis: hypoventilation	Diagnosis: Hypoventilation	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800	

Code	Description	Text message (Note: ' _ ' represents the space character 20H)	New with version	Used by device	Remarks
7FH	SmartCare diagnosis: insufficient ventilation	Diagnosis Insufficient Ventilat.	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800	
80H	SmartCare diagnosis: tachypnea	Diagnosis: Tachypnea	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800	
81H	SmartCare diagnosis: severe tachypnea	Diagnosis: Severe Tachypnea	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800	
82H	SmartCare diagnosis: unexplained hyperventilation	Diagnosis: Unexplained Hypervent	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800	
83H	SmartCare diagnosis: central hypoventilation	Diagnosis: Central Hypovent	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800	
84H	SmartCare: patient with COPD	COPD	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800	
85H	SmartCare: patient with neurologic disorder	Neurologic Disorder	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800	
86H	SmartCare: night rest requested	Night rest	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800	
87H	Sigh added to current ventilation mode	/Sigh	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800	
88H	Cuvette type: disposable	disposable cuvette	01.01	Oxylog 3000 plus, Oxylog VE300	
89H	Cuvette type: reusable	reusable cuvette	01.01	Oxylog 3000 plus, Oxylog VE300	
8AH	Hose type: disposable	disposable hose	01.01	Oxylog 3000 plus, Oxylog VE300	

Code	Description	Text message (Note: ' ' represents the space character 20H)	New with version	Used by device	Remarks
8BH	Hose type: reusable	reusable hose	01.01	Oxylog 3000 plus, Oxylog VE300	
8CH	Hose size: adult	Hose size: adult	01.01	Oxylog 3000 plus, Oxylog VE300 (added with SW 2.00 for Oxylog VE300)	
8DH	Hose size: pediatric	Hose size: pediatric	01.01	Oxylog 3000 plus, Oxylog VE300 (added with SW 2.00 for Oxylog VE300)	
8EH	Variable PS added to current ventilation mode	/VariablePS	01.00	Evita Infinity V500, Evita V600, Evita V800	
8FH	SmartCare zone of respiratory comfort: Customization active	Customizing	01.00	Evita V300, Evita Infinity V500, Evita V600, Evita V800	
90H	Overall device check result fully functional	Dev.chk.res. fully functional	01.00	Perseus A500, Atlan, Primus Family, Zeus IE (added with SW 2.02 for Zeus IE)	
91H	Overall device check result conditionally functional	Dev.chk.res. cond. functional	01.00	Perseus A500, Atlan, Primus Family, Zeus IE (added with SW 2.02 for Zeus IE)	
92H	Overall device check result non functional	Dev.chk.res. non functional	01.00	Perseus A500, Atlan, Primus Family, Zeus IE (added with SW 2.02 for Zeus IE)	

Code	Description	Text message (Note: ' _ ' represents the space character 20H)	New with version	Used by device	Remarks
9FH	Tracheal pressure calculation enabled as real-time value (independent of ATC adjunct)	Ptrach Calculation enabled	01.02	Evita V300, Evita Infinity V500, Evita V600, Evita V800, Babylog VN500, Babylog VN600, Babylog VN800 (added with SW 2.30 for Babylog VN500, added with SW 2.30 for Evita Infinity V500)	
A0H	Date and time of last self test execution (variable text transmitted with this code. Time and date are represented as specified for the Time & Date Update Response, see "Protocol Definition Dräger RS 232 MEDIBUS" for further information. Note: The month representing letters (MMM) equal German month abbreviations.)	HH:MM:SSDD-MMM-YY	01.03	Perseus A500, Atlan, Primus Family (added with SW 1.11 for Perseus A500, added with SW 4.51 for the Primus Family)	
A1H	Date and time of last leakage test execution (variable text transmitted with this code. Time and date are represented as specified for the Time & Date Update Response, see "Protocol Definition Dräger RS 232 MEDIBUS" for further information. Note: The month representing letters (MMM) equal German month abbreviations.)	HH:MM:SSDD-MMM-YY	01.03	Perseus A500, Atlan, Primus Family (added with SW 1.11 for Perseus A500, added with SW 4.51 for the Primus Family)	
A2H	Warm-up mode is active	Warm-up Mode	01.09	Babyleo TN500, Babyroo TN300	
A3H	Weaning mode is active	Weaning Mode	01.09	Babyleo TN500	
A4H	Tolerate cooling procedure is active	Cooling Mode	01.09	Babyleo TN500, Babyroo TN300	
A5H	Closed care therapy is active	Closed Care Therapy	01.09	Babyleo TN500	
A6H	Multi-step recruitment is active	Multi-step recruitment	01.10	Perseus A500, Atlan, Zeus IE (added with SW 2.00 for Atlan, added with SW 2.00 for Perseus A500, added with SW 2.02 for Zeus IE)	

Code	Description	Text message (Note: ' ' represents the space character 20H)	New with version	Used by device	Remarks
A7H	One-step recruitment is active	One-step recruitment	01.10	Perseus A500, Atlan, Zeus IE (added with SW 2.00 for Atlan, added with SW 2.00 for Perseus A500, added with SW 2.02 for Zeus IE)	
A8H	CPR ventilation is paused	CPR paused	01.12	Oxylog VE300 (added with SW 2.00 for Oxylog VE300)	
A9H	Check the calibration status of the Medibus.X device as some parameters may not be calibrated and therefore the measurements might not achieve the nominal accuracy	Check Calibration	01.12	Perseus A500, Atlan, Zeus IE (added with SW 2.00 for Atlan, added with SW 2.03 for Perseus A500, added with SW 2.02 for Zeus IE)	
AAH	Active humidification, inspiratory heated	Active humid., insp. heated	01.15	Evita V600, Evita V800, Babylog VN600, Babylog VN800 (added with SW 3.00 for Evita V800, Evita V600, Babylog VN800 and Babylog VN600)	Edition 22: Starting with SW 3.00 Evita V600, Evita V800, Babylog VN600 and Babylog VN800 use this code instead of code 76H.
ABH	Active humidification, dual heated	Active humid., dual heated	01.15	Evita V600, Evita V800, Babylog VN600, Babylog VN800 (added with SW 3.00 for Evita V800, Evita V600, Babylog VN800 and Babylog VN600)	Edition 22: Starting with SW 3.00 Evita V600, Evita V800, Babylog VN600 and Babylog VN800 use this code instead of code 77H.

Infusion pump setting responses

Code	Data description	Unit	New with version	Used by device	Remarks
00H	Infusion Rate	mg/h	01.01	Zeus IE	
01H	Infusion Rate	mg/min	01.01	Zeus IE	
02H	Infusion Rate	mL/h	01.01	Zeus IE	
03H	Infusion Rate	mL/min	01.01	Zeus IE	
06H	Infusion Rate	mg/kg/h	01.01	Zeus IE	
07H	Infusion Rate	mg/kg/min	01.01	Zeus IE	
0AH	Infusion Rate	ug/kg/h	01.01	Zeus IE	
0BH	Infusion Rate	ug/kg/min	01.01	Zeus IE	
0CH	Bolus Rate	mg/h	01.01	Zeus IE	
0DH	Bolus Rate	mg/min	01.01	Zeus IE	
0EH	Bolus Rate	mL/h	01.01	Zeus IE	
0FH	Bolus Rate	mL/min	01.01	Zeus IE	
12H	Bolus Rate	mg/kg/h	01.01	Zeus IE	
13H	Bolus Rate	mg/kg/min	01.01	Zeus IE	
16H	Bolus Rate	ug/kg/h	01.01	Zeus IE	
17H	Bolus Rate	ug/kg/min	01.01	Zeus IE	
18H	Medikation Concentration	mg/mL	01.01	Zeus IE	
19H	Medication Concentration	ug/mL	01.01	Zeus IE	
1AH	Bolus Volume	mL	01.01	Zeus IE	
1BH	Bolus Volume	mg/kg	01.01	Zeus IE	
1CH	Bolus Volume	ug/kg	01.01	Zeus IE	

Infusion pump current measured data and alarm limit response

Code	Code page	Data description	Unit	LL	HL	New with version	Used by device	Remarks
00H	1	Line pressure	mmHg		X	01.01	Zeus IE	
01H	1	Total infused volume	mL		X	01.01	Zeus IE	
02H	1	Time rest	min			01.01	Zeus IE	

Infusion pump alarm status response

Code	Code page	Alarm description	Alarm phrase	New with version	Used by device	Remarks
00H	1	base related (IV-SYSTEM COM FAILURE or IV-SYSTEM INCOMPATIBLE)	IV SYS ERROR	01.01	Zeus IE	
01H	1	pump related (PUMP REMOVED)	IV REMOVED	01.01	Zeus IE	
02H	1	pump related (INFUSION X STOPPED/SYRINGE X UNLOCKED)	IV INOP	01.01	Zeus IE	
03H	1	pump related (INF. LINE X BLOCKED)	IV OCCL	01.01	Zeus IE	
04H	1	pump related (PUMP X INFUSION LIMIT, the volume to infuse has reached)	IV HIGH	01.01	Zeus IE	
05H	1	pump related (PUMP X BATTERY?/PUMP X FAILURE)	IV ERR	01.01	Zeus IE	
06H	1	pump related (IV System Safety/PUMP X UNLOCKED)	IV COM?	01.01	Zeus IE	
07H	1	pump related (PUMP X EMPTY SYRINGE MODE)	SYR EMPTY	01.01	Zeus IE	

Infusion pump text messages

Code	Description	Text message	New with version	Used by device	Remarks
00H	Power on	Power on	01.01	Zeus IE	
01H	Bolus is running	Bolus is running	01.01	Zeus IE	
02H	Infusion is running	Infusion is running	01.01	Zeus IE	
03H	Device is booting	Device is booting	01.01	Zeus IE	
04H	Pump is accepted	Pump is accepted	01.01	Zeus IE	
05H	Pump is connected	Pump is connected	01.01	Zeus IE	
06H	Pump is disabled	Pump is disabled	01.01	Zeus IE	
0BH	Pump stopped	Pump stopped	01.01	Zeus IE	
0DH	Drug short name (variable text string transmitted with this code, see "MEDIBUS.X, Rules and Standards for Implementation" for further information)		01.01	Zeus IE	
0EH	Drug long name (variable text string transmitted with this code, see "MEDIBUS.X, Rules and Standards for Implementation" for further information)		01.01	Zeus IE	
0FH	Dosing Unit (variable text string transmitted with this code, see "MEDIBUS.X, Rules and Standards for Implementation" for further information)		01.01	Zeus IE	
11H	Concentration Unit (variable text string transmitted with this code, see "MEDIBUS.X, Rules and Standards for Implementation" for further information)		01.01	Zeus IE	

Drägerwerk AG & Co. KGaA

Moislinger Allee 53 – 55
23542 Lübeck, Germany
Tel +49 451 882-0
Fax +49 451 882-2080
www.draeger.com



9052608 – PD 6500.800 en

© Drägerwerk AG & Co. KGaA

Edition: 22 – 2024-03

(Edition: 1 – 2010-06)

Dräger reserves the right to make modifications
to the equipment without prior notice.